


Dental Instruments Catalogue



LEXER
International

Manufacturers, & Exporters of Dental, Surgical,
Orthodontic, Hollow Wears & Manicure Instruments.

IMPORTANT INFORMATION**QUALITY SERVICES & CUSTOMER SATISFACTION**

To provide "Excellent services and economically prices" is the motto of our business. All instruments are made of stainless steel other wise specified.

All instrument are stamped with words "stainless Pakistan or any other country to whom manufactured buyers own brand may also be engraved on demand, "Die-Sunk" procedure is generally used.

Special items can also be supplied if sketch, drawing or samples are provided to us or refer any reference from reputed catalogue for our guidance.

All kind of Surgical instruments of stainless steel made the specific material according to its use.

Specification of material used in production.**Cutting Instruments**

For (scissors, Cutter)	Stainless steel AISI-410/420 (Rock-well hardness-48 to 58) Hot drop forged
For (Forceps, Holding Clamp)	Stainless steel AISI-410/420 (Rock-well hardness-42 to 48) Hot drop forged
For (Scalpel, Handles)	Stainless steel AISI-420 (Rock-well hardness-48 to 58) Hot drop forged
For (Speculums, Retractors)	Stainless steel Sheet-300 Series 18/8 (18% Chromium & 8% Nickel)

Please Note:

The delicate, hand-crafted surgical instrument use only for its intended purpose to Ensure effectiveness and maximum life. Proper care must be given during cleaning, sharpening and sterilization to preserve design and integrity and to prolong usefulness. Do not expose instruments to sterilization temperatures above 350° F/177° C or to higher than recommended chemical concentrations." Longer Life for Surgical Instrument" for proper cleaning, sterilization and other instrument care.

Inspected the instrument after each use for: pitting; corrosion; discoloration; broken, loose or bent components; cracked handles; dull cutting edges; over sharpened or worn out blades; chipped or peeled plating; damaged hinges, beaks or blades. Instruments found to Have these or other compromising conditions should be removed from service.


QUALITY MANAGEMENT SYSTEM

LEXER INTERNATIONAL. is an employee owned company of Craftsmen and has got a good reputation for quality and reliability. Our significant manufacturing steps are being described here-under briefly:


1. **MATERIAL TESTING:** The purchased raw material is thoroughly tested under PQA. authority. A rigid quality control policy is being adopted before releasing to production. This test is taken to check the standard of steel as per desired specifications and is considered the most important factor of our business.
2. **FORGING & HAMMERING:** After testing, the stainless steel strips/bars are divided into different sizes as per our requirements and our highly experienced craftsmen are specialized in forging (Hot/Cold Drop) to perform all standard press operations in forming and deep drawing of all type of instruments. We anneal the forging and prepare for the milling process.
3. **MACHINE & TOOL WORK:** Our highly skilled technicians have a great knowledge about machine milling process. We form accurately male/female portion of box-joint forceps, as well as, serrations/ratches in uniformity. We are equipped to mill, drill, bore turn and grind all type of metals. We have in house facility of custom tool and die production, maintenance and repair.
4. **FILLING:** We are specialized in handmade fabrication. First, we assemble male/female portion after setting/adjusting of Extracting Forceps and Scissors, tips, serrations, ratches in uniformity and file precisely to make it semi-finish under strictly quality control specifications.
We are now moving from basic hand manufacturing techniques to mechanization to production run.
5. **HEAT TREATMENT:** All instruments are being tempered maintaining "Rock-Well Hardness" for Extracting Forceps & Retractor 45° to 50° degree and Scalar, Gauges and Scalpel 50° to 55° degree. A test is taken to check the hardness after tempering of each lot to maintain the uniformity and, before proceeding further, we clean all instruments by various chemicals. In addition, we will buy specialized Q/C equipment in near future to meet special customer Q/C. Specifications.
6. **RIVETING & WELDING:** We perform a wide variety of assembly operations from prototype quantities to production run. We rivet precisely each instrument such as, Extracting Forceps & Retractor, etc. But, for some kind of instruments, we use spot welding.
7. **GRINDING & POLISHING:** For precise grinding and excellent finish, we specialize in applying polishing wheels to belts to improve the overall finish in high uniformity utilizing best of polishing luster. We use electro-plating plants to polish serrations/ratches and joints of stainless steel instruments. For satin/matt finish instruments, we use glass heads to achieve uniformity. All instruments are being polished under rigid quality control specifications.
8. **PASSIVATION:** This method actually provides the instruments with its corrosion resistance properties. Passivation, which follows the final polishing step, is a chemical bath, which creates an oxidized layer on the surface of instruments. We are now equipped with most modern process.
9. **ULTRASONIC CLEANING:** We are full equipped with most modern ultrasonic cleaning plants imported from "Branson Cleaning Equipment Co. USA", and all instruments are being cleaned properly utilizing various chemicals such as, Tri-Chloroethylidine etc.
10. **INSPECTION & STAMPING:** We inspect all instruments under rigid Q/C. specifications and give the assurance to our relative clients that all instruments shall be free from all defects such as, burrs, sharp edges, cracks, pitting, scratches, grease / excessive oil and foreign matters etc. We also guarantee for its workmanship. For stamping, we use both process i.e. die sunk and electrolytic marking and are being marked under customer's specifications.
11. **PACKING, LABELING AND BAR CODING:** We pack all instruments under strictly Q/C. specifications as per customer's demand. Normally, we pack individually in poly bag and further 12 pieces in one box mentioning outside with Part Number, Description and Lot Number. All instruments are packed in poly bags/plastic/colour bags, Blister Packing with adequate desiccant in each poly bag to absorb moisture. We are specialized in packaging, labeling and bar coding as per customer demand.

QUALITY MANAGEMENT SYSTEM

Process Model



PROCESS FLOW CHART






01 EXTRACTING FORCEPS	01 To 16	
02 ROOT ELEVATORS	17 To 38	
03 DIAGNOSTIC, PERIODONTAL & RESTORATIVE INSTRUMENTS	39 To 92	
04 MAXILLOFACIAL SURGERY INSTRUMENTS	93 To 127	
05 ORTHODONTIC PLIERS & INSTRUMENTS	128 To 139	
06 DENTAL IMPLANT INSTRUMENTS	140 To 170	
07 IMPRESSION TRAYS & HOLLOWWARE	171 To 179	
08 INDEX	180 To 183	

INDEX



EXTRACTING FORCEPS

Dear Valued Buyers if any Universal Fig. No. of Extracting Forceps not Found in our Catalogue you Can Refer the Fig. No. From These Diagram.



Dear Valued Buyers these are our Handle Styles



P
PREMIUM



E
ERGONOMIC



S
STANDARD

What Is The Function of Extracting Forceps : To Remove Teeth From Bony Socket**Extraction Forceps**

Dental forceps are pliers-like instruments that are used to grasp teeth.

A dentist will usually have a number of different ones on hand, each having a design that's tailored to:

The specific shape of the tooth it's intended to remove (such as incisors vs. molars having two roots vs. molars having three roots).

The location of the tooth in the mouth (such as front vs. back, or left lower molar vs. right upper molar.)

How they're used.

A dentist will grasp a tooth with their forceps and then firmly and deliberately rock it back and forth as much as it will.

Because the bone that surrounds the tooth is compressible, the socket will gradually start to expand.

As it does, the range of the dentist's side-to-side motions will increase.



In addition to this rocking motion, the dentist will also rotate the tooth back and forth.

This twisting action helps to rip and tear the tooth away from the ligament that binds it in place.



At some point, the socket will be enlarged enough, and the ligament torn away enough, that the tooth can be easily removed.



EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملفظ (نمط اللغة الإنجليزية)

Used to remove teeth from alveolar bone







L1-001
1
Upper Anteriors



L1-002
7
Upper Premolars



Straight extracting forceps for upper incisors and canines



L1-003
17
Upper Molars Right






L1-004
18
Upper Molars Left



L1-005
18A
Upper Molars



L1-006
29
Upper Roots



L1-007
29N
Upper Roots Narrow



L1-008
30
Upper Roots






L1-009
55
Separating Upper Molars

EXTRACTING FORCEPS (English Pattern)



- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forma Inglesa)
- استخراج الملفظ (نمط اللغة الإنجليزية)



L1-010
41
Upper Roots





L1-011
44
Upper Roots



L1-012
51
Upper Roots



L1-013
51A
Upper Roots



L1-014
67A
Upper Third Molars



L1-015
76N
Upper Roots



L1-016
89
Upper Molars Right



L1-017
90
Upper Molars Left



L1-018
56
Separating Lower
Molars





Extracting forceps for separating upper last molar



EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملفظ (نمط اللغة الإنجليزية)






L1-019
13
Lower Premolars



L1-020
22
Lower Molars



Curved extracting for lower premolars






L1-021
33
Lower Roots



L1-022
33A
Lower Roots



L1-023
33L
Lower Roots



L1-024
73
Lower Molars



L1-025
74
Lower Anteriors
and Roots



L1-026
74N
Lower Roots



L1-027
75
Lower Premolars

EXTRACTING FORCEPS (English Pattern)



- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)

- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملفظ (نمط اللغة الإنجليزية)

A-TRAUMATIC EXTRACTION

Improved tapered beaks for complete root contact. Cross serrations.
Deep Gripe Extracting Forceps



• No Fracture

• No Tissue Damage

• Extraction With Easy



Thin Tip

L1-028
34Upper incisors and
premolars
deep-gripping

Wide Tip

L1-029
34AUpper centrals / canines
deep-gripping

Thin Tip

L1-030
Fig. 35upper permolars
deep-gripping

Wide Tip

L1-031
35AUpper incisors and
premolars
deep-gripping

Thin Tip

L1-032
36Lower premolars
deep-gripping

Wide Tip

L1-033
36ALower premolars
deep-gripping**L1-034**lower permolars
deep-gripping**L1-035**
F-7C

lower molar

**L1-036**
9Apex Forceps Lower Universal
Deep Root Pick**L1-037**
79NLower third molars
deep-gripping**L1-038**
79

Lower Third Molars

Ogden-Felsch extracting forceps
for upper third molars




EXTRACTING FORCEPS (English Pattern)


- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملفظ (نمط اللغة الإنجليزية)




L1-039
31
Lower Roots




L1-040
79
Lower Third Molars




L1-041
45
Lower Roots




L1-042
86
Lower Molars





L1-043
86A
Lower Molars



L1-044
86B
Lower Molars




L1-045
86C
Lower Molars



Extracting forceps for lower molars



L1-046
87
Lower Molars



EXTRACTING FORCEPS (English Pattern)



- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملفظ (نمط اللغة الإنجليزية)



EXTRACTION FORCEPS # ATRAUMATIC

Forceps "Blade Beaks"


1. "Blade Beaks" smoothly penetrate in contact with the root.
2. They have large contact surface for a secure grip.
3. They take care of the tissues and assure tooth fracture risk reduction.
4. They give better control over the root during the extraction operations.



L1-047
1
Upper Centrals & Canines




L1-048
7
Upper Premolars



L1-049
13
Lower Premolars



L1-050
51A
Upper Roots



L1-051
74N
Lower Roots



Cryer extracting forceps for upper incisors, premolars and roots



L1-052
79
Lower Third Molars



EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملفظ (نمط اللغة الإنجليزية)

EXTRACTION FORCEPS # SECURE

Secure-Line Extraction Forceps

These forceps will allow you to complete any extraction and work more efficiently even if complications arise. If teeth break during extraction, these are the tools you need to rapidly, safely solve your problems. These forceps have special, all-purpose, jaws designed to cover a broad range of applications and have sharp pointed tips that allow separating the tips of roots from surrounding tissues by lightly twisting them and then securely trapping the extracted roots in the gaps between their jaws.



Extracting forceps for upper roots



L1-053
Lower Roots




L1-054
Lower Roots




L1-055
Lower Roots




L1-056
Lower Roots



L1-057
Upper Roots



L1-058
Upper Roots






L1-059
Upper Roots

EXTRACTING FORCEPS (English Pattern)



- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forma Inglesa)
- استخراج المفتق (نمط اللغة الإنجليزية)



L1-060
Upper Roots



L1-061
Upper Roots



L1-062
Upper Roots



L1-063
Upper Roots



Straight extracting forceps for upper roots



L1-064
300
Upper roots
Serrated



L1-065
301
Lower roots
Serrated



L1-066
33
Lower Roots
European Style
Serrated



L1-067
51
Upper Roots
European Style
Serrated



L1-068
97
Upper Roots
European Style
Serrated



L1-069
233
Lower Roots
European Style
Serrated



L1-070
74
Lower Roots
European Style
Serrated



L1-071
74N
Lower Roots
European Style
Narrow Beaks
Serrated



EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملفظ (نمط اللغة الإنجليزية)

TOOTH FORCEPS - WISDOM



L1-073
551
Lower Molars and
Wisdoms, Right



L1-074
552
Lower Molar and
Wisdom Roots, Left



L1-075
552
Lower Molar and
Wisdom roots, Right

ANATOMIC HANDLE

L1-072
550
Lower Molars and
Wisdoms, Left



L1-076
A-51L
Upper Roots

L1-077
46
Lower Roots



L1-078
97
Upper roots & fragments



Routurier extracting forceps for lower molars
and third molars, right




Angled extracting forceps for upper roots

EXTRACTING FORCEPS (English Pattern)



- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملفظ (نمط اللغة الإنجليزية)

TOOTH FORCEPS





EXTRACTING FORCEPS (English Pattern)

LEXER
International

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملفظ (نمط اللغة الإنجليزية)



L1-090
English Pattern
Klein Fig. 4
Pedodontic Upper root



L1-094
English Pattern Klein Fig. 5
Pedodontic Lower Incisors.



L1-091
English Pattern
Klein Fig. 2
Pedodontic Upper Premolars




L1-093
English Pattern
Klein Fig. 1
Pedodontic, Upper Incisors



L1-095
English Pattern
Klein Fig. 6
Pedodontic Lower Molars



L1-092
English Pattern
Klein Fig. 3
Pedodontic Upper Molars Universal



Baby Tooth Extraction Forceps

L1-097
Set of 7 Pieces



L1-096
English Pattern
Klein Fig. 7
Pedodontic Lower Molars,

EXTRACTING FORCEPS (American Pattern)



- Щипцы экстракционные (американская форма)
- Daviers - forme américane
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzepes per estrazione - forma americana
- سخراج الملقط (النمط الأمريكي)



L1-098
Fig. W1
Winter
upper centrals
and canines



L1-099
Fig. 1
upper incisors
and roots



L1-100
Fig. 2
upper and
lower roots



L1-101
4
upper and lower premolars



L1-102
Fig. 10 H
upper molars



extracting forceps for upper left molars



EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Daviers - forme américane
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzepes per estrazione - forma americana
- سخراج الملقط (النمط الأمريكي)



L1-103
Fig. 15
Harris
lower molars



L1-104
Fig. 16
lower molars




extracting forceps for upper incisors,
premolars and roots



L1-105
Fig. 17
lower molars



L1-106
Fig. 18R
Harris
upper molars, right



L1-107
Fig. 23
lower molars



L1-108
Fig. 32 Parmly
upper premolars
and molars



L1-109
Fig. 18 L
Harris
upper molars, left

EXTRACTING FORCEPS (American Pattern)



- Щипцы экстракционные (американская форма)
- Daviers - forme américane
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzepor per estrazione - forma americana
- سخراج الملقط (النمط الأمريكي)



L1-110
53 L
upper molars, left



L1-111
53 R
upper molars,
right



L1-112
65
upper incisors and roots



L1-113
69 Tomas
upper and lower roots



L1-114
88L Nevius
upper molars, left



L1-115
88R Nevius
upper molars, right



extracting forceps for
upper third molars



L1-116
150
upper incisors,
premolars, roots



L1-117
151
Lower incisors,
premolars, roots



L1-119
151 S
lower teeth,
for children



L1-121
222
Lower third molars

L1-120
210
upper third molars



EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Daviers - forme américane
- Fórceps de extracción - forma americana
- Щипці екстракційні (американська форма)
- Pinzepes per estrazione - forma americana
- سخراج المقلظ (النمط الأمريكي)



Curved extracting forceps for fractured or decayed lower teeth

L1-123
Lower Molar
Fig. 23 Cow Horn



L1-124
Fig. 74
Lower Anterior



L1-125
Fig. 88L
Upper Molar Left



L1-126
Fig. 88R
Upper Molar Right



L1-127
Fig. 74
Lower Anterior



L1-128
Fig. 106
Lower Incisors



L1-129
Fig. 174
Extra Narrow
Lower Anterior
Apical Retention
Forceps



L1-122
Fig. 1
Upper Anterior



L1-130
Fig. 859
Lower Molars
Roots



DENTAL SYRINGES



- Шприцы стоматологические карпульные
- Seringues à usage dentaire
- Dentalspritzen
- Шприці стоматологічні карпульні
- Jeringas dentales
- المحقن الأنفان



L1-131 1.8ML
L1-132 2.2ML



L1-133 1.8ML
L1-134 2.2ML



L1-135 1.8ML
L1-136 2.2ML



L1-137 1.8ML
L1-138 2.2ML



L1-139 1.8ML
L1-140 2.2ML



L1-141 1.8ML
L1-142 2.2ML



L1-143
Dental Capsule Applier Applicator
Gun for (GC Fuji) (SDI Rive)



L1-144
1001
Intraligamental
Anesthesia Syringe
in Black Colour



L1-145 1.8ML
L1-146 2.2ML
Angled



L1-147
Aluminium - Dark 1.8 ml



L1-148
Aluminium - Silver 1.8 ml



L1-149
Aluminium - Orange 1.8 ml



L1-150
Aluminium - Green 1.8 ml




L1-151
Aluminium - Blue 1.8 ml



ROOT ELEVATOR

LEXER
International



HANDLES FOR ROOT ELEVATORS POINTS A, B, C, D, E, F



PREMIUM QUALITY DENTAL INSTRUMENTS



Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.

Dental Root Elevators

To look at one of these instruments you might mistake it for a narrow screwdriver.

That's because just like one an elevator has a handle and then a specially designed "blade" or tip portion.

How they're used.

The tip of an elevator is designed to be wedged into the ligament space between a tooth and its surrounding bone (as shown in our animation).

As the elevator is forced into and twisted around in this space, the tooth is rocked around and pressed against the walls of its socket.

This helps to both expand the shape of the socket and separate the tooth from its ligament.

As this work is continued, the tooth gradually becomes more and more mobile.

In some cases, the dentist may be able to completely remove the tooth just using their elevator.

If not, they will choose a point where they feel they have accomplished as much as they can with it and will then switch to using extraction forceps.





ROOT ELEVATOR

- Элеватор корневой экстракционный
- Élevateurs à racines
- Wurzelheber
- Слеватор кореневий экстракційний
- Elevadores para raíces
- المصاعد الجذر

Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



L2-001
1
Apical 301



ROOT ELEVATOR



- Элеватор корневой экстракционный
- Élevateurs à racines
- Elevadores para raíces
- المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



L2-016
301
Apexo

For displacement of deeply seated roots



L2-017
303
Apexo



L2-018
301W
Seldin



L2-019
304W
Seldin



L2-020
1
Heidbrink
Straight



L2-021
2
Heidbrink
Straight

Upper and lower anterior teeth and roots



upper root elevator



upper root elevator



L2-022
80



L2-023
81



L2-024
81A



L2-025
81E



L2-026
40
Heidbrink
Straight



L2-027
41
Heidbrink
Straight

Upper anterior roots and small teeth For displacement of anterior roots

Upper and lower anterior teeth and roots




ROOT ELEVATOR

- Элеватор корневой экстракционный
- Élevateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر



Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



L2-038
Straight
3mm

ROOT ELEVATOR




- Элеватор корневой экстракционный
- Élevateurs à racines
- Wurzelheber

- Сліватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.





ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



CRYER

Cryer root or fragment elevator



ROOT ELEVATOR



- Элеватор корневой экстракционный
- Élevateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



SELDIN

Seldin root or fragment elevator



Upper and lower molar roots





ROOT ELEVATOR

- Элеватор корневой экстракционный
- Élevateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



L2-096
4
Schmeckebier
Apexo

Schmeckebier and Miller designs used for displacement of broken or deep-seated root tips.



L2-097
5
Schmeckebier
Apexo



L2-098
71
Miller
Apexo



L2-099
72
Miller
Apexo



L2-100
6
Potts

For impacted upper 3rd Molars without bony obstruction.



L2-101
7
Potts



L2-102
73
Miller
Apexo



L2-103
74
Miller
Apexo



L2-104
19



L2-105
20

Shorter leverage for lower 1st and 2nd Molar roots. Also useful for premolar and molar deciduous teeth.



lower roots elevator

ROOT ELEVATOR



- Элеватор корневой экстракционный
- Élevateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



L2-106
38
Serrated

For removal of root fragments

L2-107
39
Serrated

L2-108
190

L2-109
191

L2-110
15
Woodward
Lower Molar roots

L2-111
16
Woodward
Lower Molar roots



ROOT ELEVATOR

- Элеватор корневой экстракционный
- Élevateurs à racines
- Elevadores para raíces
- المصاعد الجذر

Function of Root Elevators:

To loosen tooth or root from bony socket prior to placement of the extraction forceps.



L2-112
23



Upper root elevator with helicoidal tip



L2-113
H5
Hourigan
Serrated



L2-114
H6
Hourigan
Serrated



L2-115
1C
Coupland
Gouge



L2-116
2C
Coupland
Gouge
Lower Molars



L2-117
3C
Coupland
Gouge
Lower 3rd Molars



L2-118
69W
Woodward Hand
Pressure Chisel



L2-119
31F
Friedman
Craig
Gouge Serrated



L2-120
32F
Friedman
Craig
Gouge Serrated



L2-121
11M
MacMillan
Gouge



L2-122
12M
MacMillan
Gouge



L2-123
13M
MacMillan
Gouge

For Upper and Lower Cuspids

ENGLISH PATTERNS



L2-124
Uppe and Lower
3rd Molars



L2-125
A1C
Coupland
Gouge



L2-126
A2C
Coupland
Gouge



L2-127
A3C
Coupland
Gouge



L2-128
A30
Cryer



L2-129
A31
Cryer

ROOT ELEVATOR



- Элеватор корневой экстракционный
- Élevateurs à racines
- Elevadores para raíces
- المصاعد الجذر

ERGO DYNAMIC ELEVATORS

Ergo Dynamic Elevators offer an alternative to standard size elevator handles.



L2-130
301
Apexo

L2-131
301
Apexo
Ergo Dynamic

L2-132
#77R Elevator
Serrated
Ergo Dynamic

L2-133
#B Cogswell
Elevator
Ergo Dynamic

L2-134
#34S
Seldin Elevator
Straight
Ergo Dynamic

L2-135
Luxating Elevator
Straight
3mm
Ergo Dynamic

L2-136
Luxating Elevator
Curved
3mm
Ergo Dynamic

L2-137
Luxating Elevator
Straight
5mm
Ergo Dynamic

L2-138
Luxating Elevator
Curved
5mm
Ergo Dynamic





ROOT ELEVATOR

- Элеватор корневой экстракционный
- Élevateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

Laxutor Root Elevators Set of 9 Pcs



ROOT ELEVATOR



- Элеватор корневой экстракционный
- Élevateurs à racines
- Elevadores para raíces
- المصاعد الجذر

Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



Heidbrink universal mini-elevator for fragments

HOWARD

Slightly heavier design than Heidbrink picks.

HEIDBRINK

APICAL



ROOT ELEVATOR

- Элеватор корневой экстракционный
- Élevateurs à racines
- Wurzelheber
- Сліватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر



L2-159
Heidbrink
Left



L2-158
Heidbrink
Straight



L2-160
Heidbrink
Right



L2-162
Heidbrink
Left



L2-163
Heidbrink
Right



L2-161
Heidbrink
Straight



L2-164
Warwick-
James
Straight



L2-165
Warwick-
James
Left



L2-166
Warwick-
James
Right



L2-167
Fig. 1
190/15



L2-168
Fig. 2
190/16



L2-169
Fig. 3
190/17



universal mini-elevator for fragments

Coupland



L2-170
192/1
Left



L2-171
192/2
Right



L2-172
192/3
Str.



L2-173
Fig. 1



L2-174
Fig. 2



L2-175
Fig. 3

ROOT ELEVATOR



- Элеватор корневой экстракционный
- Élevateurs à racines
- Elevadores para raíces
- المصاعد الجذر



L2-176
Fig. 1
Chompret



L2-177
Fig. 2
Chompret



L2-178
Fig. 3
Chompret



L2-179
Fig. 4
Chompret



Lecluse elevator for lower third molar dislocation

CROSS-BAR HANDLE



L2-180
11LX

L2-181
11RX



Chompret syndesmotome, sickle pattern



L2-182
12LX



L2-183
12RX



L2-184
14LX



L2-185
14RX



L2-187
1RX



ROOT ELEVATOR

- Элеватор корневой экстракционный
- Élevateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

Function of Root Elevators: To loosen tooth or root from bony socket prior to placement of the extraction forceps.



L2-188
775/3
Winter
11 L



L2-189
775/4
Winter
11 R



L2-190
775/5
Winter
12 L



L2-191
775/6
Winter
12 R



L2-192
775/7
Winter
13 L



L2-193
775/8
Winter
13 R



L2-194
winter
Left

L2-195
winter
right



L2-196
775/1
Winter
1 L

L2-197
775/2
Winter
1 R

L2-198
Winter
1L


L2-199
Winter
1R

ROOT ELEVATOR

- Элеватор корневой экстракционный
- Élevateurs à racines
- Wurzelheber

- Сліватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

LUXACUT



Luxacut Luxation Instrument

Over recent years, there has been an increased demand for reduced trauma during extraction procedures. The use of elevators has become widespread for the loosening of teeth before the application of forceps. Luxation instruments tips are more finely tapered than elevators tips. Their characteristics are a combination of a Periotome and an elevator. The sharp tip enables more efficient cutting of the periodontal ligament, whilst the tapering shaft assists with controlled expansion of the bony socket. The extraction operation may be performed with the minimum amount of tissue damage. Healing is therefore quicker and the operation is less distressing for the patient and surgical team alike.


The Luxacut Kit comprises of seven Luxation instruments, a straight, curved, inverted curved & contra blade version in Three different widths. Also available in dual edge blades

- **L2-196-L2S** Straight Blade- Green-2 mm- Apical and/or Interproximal
- **L2-198-L3C** Curved Blade- Black- 3 mm General use/Starter
- **L2-197-L3S** Straight Blade- Grey- 3 mm Interproximal
- **L2-200-L3CA** Contra Angle- Blue- 3 mm Lingual and/or distal-Molars
- **L2-199-L3IC** Inverted Curved- Light Blue- 3 mm- General, Lingual and/or Distal
- **L2-202-L5C** Curved Blade- Brown- 5 mm- Large molar roots-General use
- **L2-201-L5S** Straight Blade- Yellow- 5 mm- Interproximal
- **L2-203-DE3** dual edge/straight grey 3 mm/1.5 mm
- **L2-204-DE5** dualedge/straight brown 5 Deep fractured/Decayed Roots

Luxacut Kit



**L2-205
LWC7**



**L2-206
LGWC7**



ROOT ELEVATOR

- Элеватор корневой экстракционный
- Elévateurs à racines
- Elevadores para raíces
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

L2-207

LUXATING ELEVATORS SET OF 10 PCS



L2-218

LUXATING ELEVATORS SET OF 10 PCS - BLUE TITANIUM COATED TIPS



ROOT ELEVATOR



- Элеватор корневой экстракционный
- Élevateurs à racines
- Elevadores para raíces
- المصاعد الجذر


L2-229

LUXATING ELEVATORS SET OF 10 PCS - GOLD COATED TIPS



L2-240

LUXATING ELEVATORS SET OF 08 PCS - BLUE TITANIUM COATED TIPS





ROOT ELEVATOR

- Элеватор корневой экстракционный
- Élevateurs à racines
- Elevadores para raíces
- المصاعد الجذر

L2-249

PERIOTOME LUXATING ELEVATORS SET OF 04 PCS - GOLD COATED TIPS



L2-250
PERIOTOME GOLD TIP - L3S



L2-251
PERIOTOME GOLD TIP - De3



L2-252
PERIOTOME GOLD TIP - L2S



L2-253
PERIOTOME GOLD TIP - L3CA



L2-254 (Double Curved 1.5mm)
L2-255 (Double Curved 2mm)
L2-256 (Double Curved 3mm)
L2-257 (Double Curved 4mm)
L2-258 (Double Curved 5mm)



L2-259
(Inverted Curved 3mm)



L2-260
(Dual Edge 3mm / 1.5mm)





L2-261
(Dual Edge 5mm / 3mm)

ROOT ELEVATOR



- Элеватор корневой экстракционный
- Élevateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر

PDL ELEVATORS





ROOT ELEVATOR

- Элеватор корневой экстракционный
- Élevateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raíces
- المصاعد الجذر



Area: distal










Area: mesial



L2-264

PDL Elevators Set of 10 Pcs



L2-265

Luxa-Tool, straight, 3.0 mm,

L2-266

Luxa-Tool, curved, 3.0 mm,

L2-267

Proximator, curved, 2.5 mm,
distal elevator,

L2-268

Proximator, curved, 2.5 mm,
mesial elevator,

L2-269

Proximator, curved, 4.5 mm,
distal elevator,

L2-270

Proximator, curved, 4.5 mm,
mesial elevator,

L2-271

Syndesmotome, curved, 4.0 mm,
for distal access

L2-272

Syndesmotome, curved, 4.0 mm,
for mesial access

L2-273

Syndesmotome, curved left,
4.0 mm, for lingual / buccal luxation,

L2-274

Syndesmotome, curved right,
4.0 mm, for lingual / buccal luxation,

ROOT ELEVATOR



- Элеватор корневой экстракционный
- Élevateurs à racines
- Wurzelheber
- Слеватор кореневий экстракційний
- Elevadores para raíces
- المصاعد الجذر

PDL ELEVATORS

Delicate Arrow tip

4,0mm
straight
curvedLuxivators straight
atraumatic luxation of the tooth2,5mm
straight
3,5mm
straight
4,5mm
straightLuxivators curved
atraumatic luxation of the tooth2,5mm
curved
3,5mm
curved
4,5mm
curved

L2-275
PDL Elevators
Set of 8 Pcs



LEXER Dental Hand Instruments

The types of hand instruments can be categorized based on their function:

- Examining the mouth and teeth
- Scaling or professional cleaning of teeth
- Cutting teeth and removing cavities
- Placing and condensing filling materials
- Carving and finishing dental fillings
- Miscellaneous types
- Instruments used for examining the mouth and teeth

Mouth mirrors are used to examine the condition of the teeth and mouth, as well as to retract the tongue and cheeks upon examination.

Probes— There are many probes available but the main types are the straight probe, Briault probe and periodontal probe.

The straight probe is sharp and it is used to examine cavities in the dentine region as well as to check the margins of restorations or fillings.

The Briault probe is also sharp and it is used to detect the junction between the enamel and dentine. Briault can also be used to detect tartar hidden in the gum pockets.

A periodontal probe can be used to measure the depth of the gum pockets and for other measurements. Instruments for cutting and scraping tissues and deposits

Scalers remove tartar (calculus) on the tooth surface as well as those hidden in the gum pockets. They can also be used to remove other surface deposits and temporary crowns.

Excavators are used to remove softened dentine and temporary fillings. The back of the blade can also be used for placing linings between filling materials. Sometimes excavators can be used for carving metal or amalgam fillings.

Chisels are used to remove unsupported enamel during cavity preparation. Unsupported enamel can be harmful for it can compromise the longevity of the filling if not removed.

Hatchets and hoes have similar purposes to chisels. Instruments used for placing and condensing filling materials Plastic instruments convey and shape filling materials and are usually flat blades that do not use heavy pressure.

Some instruments have rounded ends for shaping and burnishing the surfaces of fillings. Plastic instruments are usually made of stainless steel for most materials though there are also thin Teflon-coated or titanium nitride instruments which are non-stick for composite fillings. **Condensers or pluggers** are used for compressing and forming filling materials. Heavy pressure is applied on use therefore these instruments are usually used for metal or amalgam fillings. Instruments used for carving and finishing fillings Carving and finishing instruments have sharp or semi-sharp blades of various shapes. They are used to carve materials by cutting and scraping.

VALUED BUYERS

Compare our scalers and curettes with the best in the industry and you will immediately recognize we are equal in every way but one... price. Note the deep precision knurling for better tactile feel. The same high-grade stainless steel tips, the same light weight (under 20 grams) for less hand fatigue, the same professional craftsmanship, and surpassing the highest industry standards. of course, any manufacturer can make these claims, but we guarantee and certify ours. We regularly manufacturing below articles in showing handles & pattern for better grip & perfection. We are also fully equipped & capable to develop these instruments as per your own specific requirements & demands in materials, sizes, shapes and patterns.



Symptoms of Gum and Dental Problems

Gum Problem Symptoms

- Sensitive or soft gums
- Swollen gums
- Gum recession
- Dry or shiny gums
- Change in gum color (from pink to pale or pink to red)
- Pain or burning in the gums
- Gums bleed more easily than usual
- Pus between the teeth and gums

Unhealthy



Dental Problem Symptoms

- Dry mouth
- Burning tongue
- Tooth decay
- Cavities
- Jaw problems
- Altered taste perception
- Bad breath

Healthy teeth and gums





Gingivitis





Periodontitis



Advanced periodontitis







HANDLES FOR MOUTH MIRRORS AND MIRRORS

- Ручки для зеркал
- Manches pour miroirs à bouche
- Mundspiegelgriffe
- Ручки для дзеркал
- Mangos para espejos de boca
- مقابض لمرايا الفم



Rhodium



L3-001
Fig. 3
Ø20mm



L3-002
Fig. 4
Ø22mm



L3-003
Fig. 5
Ø24 mm



Cone-Socket Rhodium



L3-004
Fig. 3
Ø20mm



L3-005
Fig. 4
Ø22mm



L3-006
Fig. 5
Ø24mm



L3-007
Mirror Handle



L3-008
Mirror Handle



L3-009
Hollow Mirror
Handle



L3-010
Hollow



L3-011
Hollow



L3-012
Mirror Handle
Stem with
Rule 35mm



L3-013
Cone-socket
Solid



L3-014
Cone-socket
Solid



L3-015
Cone-socket
Hollow

MOUTH MIRRORS



- зеркал
- Espejos de la boca
- Mundspiegel
- дзеркал
- Miroirs à la bouche
- مرآة الفم



Plain



Cone-Socket Plain

L3-016
Fig. 3
Ø20mmL3-017
Fig. 4
Ø22mmL3-018
Fig. 5
Ø24mmL3-019
Fig. 6
Ø26mmL3-020
Fig. 3
Ø20mmL3-021
Fig. 4
Ø22mmL3-022
Fig. 5
Ø24mm

Magnifying

Cone-Socket Magnifying

L3-023
Fig. 3
Ø20 mmL3-024
Fig. 4
Ø22 mmL3-025
Fig. 5
Ø24 mmL3-026
Fig. 3
Ø20mmL3-027
Fig. 4
Ø22mmL3-028
Fig. 5
Ø24mm



EXAMINATION KIT

LEXER
International

L3-029

Examination kit Set of 3 Pieces



Probe

Mirror

Meriam
16cm



L3-030

Examination Set of 4 Pieces



Exp,23/17-A

Exavators

Meriam
16cm

Mouth Mirror

L3-031

Initial examination Set of 3 Pieces



L3-032

Initial examination Set of 3 Pieces



L3-033

Initial examination Set of 3 Pieces



TWEEZERS



- Пинцеты
- Précelles - Bruxelles
- Zahnpinzelten

- Пінцети
- Pinzas para algodón
- ملقط



L3-034
London College
15cm



L3-035
Semken Taylor
Str.



L3-036
Semken Taylor
Cvd.



L3-037
Meriam
16cm



L3-038
Semken Taylor
1x2 Straight



L3-039
Semken Taylor
1x2 Curved



L3-040
Perry
13cm



L3-041
College with Lock
15cm Cvd.



L3-042
Gerald
17.5cm



L3-043
with lock
15cm



TWEEZERS

- Пинцеты
- Précelles - Bruxelles
- Zahnpinzelten
- Пінцети
- Pinzas para algodón
- ملاقط



L3-044
Adson-TC
12cm

L3-045
Adson-TC
12cm 1x2

L3-046
Adson-TC
15cm 1x2

L3-047
Adson
12cm

L3-048
Adson
12cm 1x2

L3-049
Adson
15cm 1x2

L3-050
Adson Brown
12cm

L3-051
Crane Kaplan
15cm



L3-052
Crane Kaplan
15cm

L3-053
Crane Kaplan
15cm

L3-054
Miller
15cm Str.

L3-055
Miller
15cm Cvd.


Dressing Forceps
serrated
L3-056 10.5 cm
L3-056 11.5 cm
L3-056 13 cm
L3-056 14 cm
L3-056 16 cm
L3-056 18 cm
L3-056 20 cm
L3-056 25 cm
L3-056 30 cm

Tissue Forceps
1x2 teeth
L3-057 10.5 cm
L3-057 11.5 cm
L3-057 13 cm
L3-057 14 cm
L3-057 16 cm
L3-057 18 cm
L3-057 20 cm
L3-057 25 cm
L3-057 30 cm



- Зонды
- Sondes
- Zahnsonden
- Зонди
- Sondas
- المستكشفيون

Single Ended Probe Solid Handle 5.5mm, Length 150mm



Explorers - single-ended



*Periodontal pocket probe
n.2N with a curved tip
for easier access to
forked areas*

*Explorer for
diagnostic
fig.23*



EXPLORERS

- Зонды
- Sondes
- Zahnsonden
- Зонди
- Sondas
- المستكشفيون

EXPLORERS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

Explorers-Single-End




L3-082
23
Shepherd's
Hook
Explorer
#30, #6, #7, #8

L3-083
#5, #6
Explorer



*Periodontal pocket
probe n.2N with a
curved tip for easier
access to forked
areas*

Double-Ended Explorers




L3-084
TU 17/23
Wilkins/Tufts
Explorer
#6

L3-085
23/UNC1
university of north
carolina
Explorer #31, #6, #7, #8

L3-086
1
#31, #6

EXPLORERS



L3-087
DG 16
#31, #6, #7
Thick, #6


L3-088
DG 16/23

L3-089
DG 16/17

L3-090
1
Abou-Rass
Explorers



Explorer with very extended and curved tip, same shape of Gracey Curette fig.11/12AF for easier access to forked areas



Periodontal pocket probe with two different tips: one cone shaped end (Williams), one flat end (Goldman-Fox)

2-4-6-8-10-12

3-6-8-11

1-2-3-5-7-8-9-10

3-6-8-11

3-6-9-12

3-5-8-10



L3-092
CP-8



L3-093
CP-10



L3-094
CP-11



L3-095
CP-12



L3-096
CP-18

1-2-3-4-5-6-7-8-9-
10-11-121-2-3-4-5-6-7-8-9-
10-11-12-13-14-15

3-6-8

1-2-3-5-7-8-9-10

1-2-3-5-7-8-9-10

3.5-5.5-8.5-11.5

L3-091
CP2



L3-097
CP-12

University of
North Carolina



L3-098
CP-15

University of
North Carolina



L3-099
PQ-O
Michigan



L3-100
PQ-W
Williams



L3-101
PQ-OW
Thin Williams



L3-102
Maryland
Laser Etched



L3-103
CP-11.5B
Screening
Probe



PROBES

DOUBLE-ENDED SCREENING PROBES



NOVATECH™ PROBES

Unique right angle design for improved adaptability in posterior regions.



Conventional Probe



Novatech Probe



WILLIAM PROBE SET DOUBLE ENDED WITH HOLLOW HANDLE 9.5MM, LENGTH 170MM

Proper care must be given during cleaning, sharpening and sterilization to preserve design and integrity and to prolong usefulness. Do not expose instruments to sterilization temperatures above 350 F/177 C or to higher than recommended chemical concentrations.



Canal or endodontic explorer




Explorer fig. CP-12S OMS with colored notches and round tip used to probe depth of the sac



L3-122

PERIODONTAL SET OF 10 PIECES SOLID HANDLE 6MM LENGTH 170MM




Williams Periodontal Pocket Probe




ENDODONTICS INSTRUMENTS

ANTERIOR



Spreader for lateral condensation

ENDODONTICS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM
POSTERIOR

POCKET PROBES



Spreader for lateral condensation




Plugger for vertical condensation



SPREADERS

HEAT CARRIERS






ROOT CANAL PLUGGER / SPREADER

LEXER
International

EXPLORERS INSTRUMENTS
SOLID HANDLE 6MM, LENGTH 170MM







- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

COMPOSITE/PLASTIC FILLING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

POSTERIOR



UNIVERSAL FILLING INSTRUMENTS





FILLING INSTRUMENTS


- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

COMPOSITE/PLASTIC FILLING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

BLACK

PLUGGERS/CONDENSERS

Hollenback



Universal cylindric filling
instrument with serrated tips



Ball filling instruments


FILLING INSTRUMENTS




- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

PLUGGERS/CONDENSERS

RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM



PREMIUM QUALITY DENTAL INSTRUMENTS






FILLING INSTRUMENTS

LEXER
International


- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Інструменти для реставрації
- Instrumentos para obturaciones
- ملع السکوک

RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

Tanner Filling Instruments




Back Action Filling Instruments



Solid Handle Size 4.25"

BURNISHERS

BALL




L3-220
1T
Tanner

L3-227
18

L3-228
26/27S

L3-229
27/29

BEAVERTAIL



L3-230
2
Single-End


L3-230
2
Double-End

L3-231
21

L3-232
21B

L3-233
2
Ladmore

ACORN



LADMORE



Ladmore burnisher

FILLING INSTRUMENTS




- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك



Sugarman file




BURNISHERS

L3-235
PKT-3R
Rounded
ConeL3-236
RomerowskiL3-237
3
LadmoreL3-234
2/29Universal cylindric filling
instrument with serrated tips

PERIODONTAL FILES

L3-239
11/12
Buck
Straight FileL3-240
1S/2S
Sugarman File
Mesial/DistalL3-241
3S/4S
Sugarman File
Buccal/Lingual

THARP

L3-238
9/10
Schluger
Curved FileL3-242
18
SprengelL3-243
1
TharpL3-244
2
Tharp




FILLING INSTRUMENTS

LEXER
International

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك


HATCHETS

RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM




PREMIUM QUALITY DENTAL INSTRUMENTS

ANTERIOR



CARVERS


RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM



TANNER

WALL

WARD





ENDODONTIC INSTRUMENTS

LEXER
International



Root Canal Explorer



EXCAVATORS



L3-273
Fig. 6/1

L3-274
Fig. 6/17

L3-275
Fig. 6/23



L3-276
1
Endoglick

L3-277
31L

L3-278
32L

L3-279
33L

PREMIUM QUALITY DENTAL INSTRUMENTS

FILLING INSTRUMENTS



- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)

- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

Ladmore burnisher

LADMORE



L3-280
Burnisher
Ladmore 3

L3-281
Burnisher
28-29

L3-282
Burnisher
26-30

L3-283
Burnisher
26-27

L3-284
T. Ball Burnisher
29-33 / C-1000
29-34

L3-285
Birtles Spatula
Applicator

L3-286
Ball Burnisher
28-34



HATCHET FILLING INSTRUMENTS

LEXER
International



L3-287
Burnisher
Pf121

L3-288
Burnisher
NYU-B1

L3-289
HOW .65mm

L3-290
Enamel Hatchet
1mm

L3-291
Enamel Hatchet
1.5mm

L3-292
Margin Trimmer
1mm 80 Degree

L3-293
Margin Trimmer
1mm 95 Degree

PREMIUM QUALITY DENTAL INSTRUMENTS



L3-294
Composite 4



L3-295
Composite 6



L3-296
Composite 7



L3-297
Composite 8



L3-298
Composite 9



L3-299
Composite 12




CARVERS INSTRUMENTS

- Инструменты для реставрации
- Instruments a modeler
- Modellierinstrumente
- Інструменти для реставрації
- Instrumentos para modelar
- أدوات تقطيع

*Hollenback small
amalgam carver fig.1/2*



CARVERS



L3-300
Hollenback H1/2

L3-301
Hollenback H1/2L

L3-302
Hollenback H3

L3-303
Tharp

L3-304
Thompson 3

L3-305
Interproximal
Carver


L3-306
U of minnesota 1

EXCAVATORS



- Экскаваторы
- Excavafteurs
- Excavatoren
- Єскаватори
- Excavadores
- حفارات

SPOONS



Excavator for manual removal of decay tissue






FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

DISCOID-CLEOIDS



L3-321
Amlagam
Condenser
.7mm-1mm

L3-322
Amlagam
Condenser
1mm-1.5mm

L3-323
Amlagam
Condenser
1.5mm-2mm

L3-324
Amlagam
Condenser
Oregon 1


L3-325
Amlagam
Condenser
Oregon 2

L3-326
Amlagam
Condenser
Oregon 3

L3-327
Discoid
Cleoid
1.5mm



*Excavator for manual removal
of decay tissue*



L3-328
Discoid Cleoid
2mm

L3-329
Discoid Cleoid
3mm

L3-330
Tanner 4

L3-331
Plastic 1

L3-332
Plastic 8A

L3-333
Plastic 5-7

L3-334
Gregg 4-5

EXCAVATORS



- Экскаваторы
- Excavafteurs
- Excavatoren
- Єкскаватори
- Excavadores
- حفارات





EXCAVATORS

LEXER
International

- Экскаваторы
- Excavafteurs
- Excavatoren
- Єскаватори
- Excavadores
- حفارات

EXCAVATORS

BLADES

ENGLISH PATTERNS



PREMIUM QUALITY DENTAL INSTRUMENTS

*Angular shaped excavator with rounded tip***L3-373**

The **2 Glick Endo Excavator** is a popular instrument for easy access and removal of dentin.

L3-374

The **31W Endo Excavator** features an extended length shank for easy access and removal of dentin.

L3-375

The **31L Endo Excavator** features an extended length shank for easy access and dentin removal.

L3-376

The **32L Endo Excavator** features an extended length shank for easy access and removal of dentin.

L3-377

The **33L Endo Excavator** features an extended length shank for easy access and removal of dentin.



GINGIVAL CORD PACKERS

LEXER
International

GINGIVAL CORD PACKERS Foratraumatic and accurate cord placement.



L3-378 BN1 Thin blade and rounded contour facilitates use in both thick and thin tissues without catching or dropping cord. Bilateral notch allows placement in limited access areas.



L3-379 S6 Ideal blade thickness with angle and blade shapes similar to the IPC carver.



L3-380 CSI-1 Serrated Also available non-serrated CSI-1 (GCPCSI1NS).



L3-381 113 Serrated Also available non-serrated 113 (GCP113NS).



L3-382 Balshi Small and thin non-serrated blade shape.



L3-383 1 Yardley Round non-serrated working end.



L3-384 7 Guyer Serrated Also available non-serrated 7 Guyer (GCPG7NS).




GCPGS6 accurate cord placement with minimal tissue trauma



ENAMEL INSTRUMENTS


MARGINTRIMMERS



HATCHES



Scaler for reduced interproximal spaces





HATCHET FILLING INSTRUMENTS

LEXER
International

CHISELS

0,5 mm



1

L3-395

1,0 mm



2

L3-396

1,5 mm



3

L3-397

0,5 mm



4

L3-398

1,0 mm



5

L3-399

1,5 mm



6



L3-400

0,5 mm



7



L3-401

1,0 mm



8



L3-402

1,5 mm




9




L3-403



COMPOSITE INSTRUMENTS TITANIUM BLUE



COMPOSITE INSTRUMENTS TITANIUM GOLD



Also available in non-titanium.



Scaler for reduced interproximal spaces

CARVERS

Used to carve anatomical features and trim excess materials.



L3-416 Inter-proximal Extremely thin, flexible blade; ideal for interproximal contouring.



L3-417 Inter-proximal Off-Angle Extremely thin, flexible blade. Off-Angle provides better access to posterior areas.



L3-418 1/2-3 Hollenback Universal adaptability. Ideal for placing, carving and contouring amalgam.




L3-419 3S Hollenback Design characteristics similar to 1/2 Hollenback but with slightly larger blades.



L3-420 3 Hollenback Design characteristics similar to 1/2-3 Hollenback but with significantly larger blades.



L3-421 8 Wiland Minn.IPC1 Extremely thin curved blade; ideal for adapting to interproximal surfaces.

COMPOSITE INSTRUMENTS
RESTORATIVE TITANIUM INSTRUMENTSL3-422
#3
Extra-Flex

PTNCIGFT3

L3-423
Mini 1

PTNCIGFTM1

L3-424
Micro-Mini

PTNCIPCS

L3-425
Medium
Placing/
Condensing

PTNCIPCM

L3-426
Large
Placing/
Condensing

PTNCIPCL PTNPLG0T

L3-427
OT
Tanner

PTNPLGH3

L3-428
3
Hollenback


PTNPLG5A

L3-430
Small/Medium
Contact
Forming

PTNCFIS/M

L3-431
Medium/Large
Contact
Forming


PTNCFIM/L

COMPOSITE INSTRUMENTS
RESTORATIVE TITANIUM INSTRUMENTSL3-432
AB1
Boghosian
PTNPFIAB1L3-433
AB2
Boghosian
PTNPFIAB2L3-434
W3
PTNPFIW3L3-435
8A
PTNPFI8AL3-436
4/5
Gregg
PTNPFIG4/5L3-437
Interproximal
Carver
PTNCVIPCL3-438
Interproximal
Carver, Long
PTNCVIPCLL3-439
A6
PTNPFI6L3-440
Contouring
Instrument
PTNCFIR/L

COMPOSITE INSTRUMENTS
RESTORATIVE TITANIUM INSTRUMENTS

BURNISHERS

FREEDMAN BURNISHERS

L3-441
21B

PTNBB21B

L3-442
27/29

PTNBB27/29

L3-443
2
Ladmore

PTNBBL2

L3-444
3
Ladmore

PTNBBL3

L3-445
Freedman
"Duckhead"


PTNPCCI

L3-446
Freedman
Small Contact
Forming
PTNFCISL3-447
Freedman
Large Contact
Forming
PTNFCILL3-448
Small/
Medium Ball
Burnisher
PTNBBS/M

COMPOSITE INSTRUMENTS
RESTORATIVE TITANIUM INSTRUMENTS

a full veneer bonding

retraction cord

L3-449
#1L3-450
#2L3-451
#3
Extra-FlexL3-452
Micro-MiniL3-453
Mini 1L3-454
Mini 4
Extra-FlexL3-455
#4
Extra-FlexL3-456
#5
Flexi-ThinL3-457
#6
Flexi-ThinL3-458
Mini 3
Extra-Flex

PTNCIGFT1 PTNCIGFT2 PTNCIGFT3

PTNCIPCS PTNCIGFTM1 PTNCIGFTM4

PTNCIGFT4

PTNCIGFT5

PTNCIGFT6 PTNCIGFTM1

FILLING INSTRUMENTS



- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)

- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

ENGLISH PATTERNS

L3-459
SerratedL3-460
SerratedL3-461
Serrated

CLEOID-DISCOIDS



L3-462

L3-463

L3-464

L3-465

L3-466

L3-467



DICAL APPLICATORS




L3-468

L3-469

L3-470

L3-471

L3-472
BerlinL3-473
HeidemannL3-474
HeidemannL3-475
Heidemann




GRACEY CURETTES

LEXER
International



Scaler for manual removal
of over-gingival backside
tartar

GRACEY SET OF 7 PCS HOLLOW HANDLE 9.5MM, LENGTH 170MM WITH SILICON GRIP



L3-476
GRACEY
1/2

L3-477
GRACEY
3/4


L3-478
GRACEY
5/6

L3-479
GRACEY
7/8

L3-480
GRACEY
9/10

L3-481
GRACEY
11/12

L3-482
GRACEY
13/14



JACQUETTE SCALERS

L3-483
30/33 JacquetteL3-484
31/32 JacquetteL3-485
34/35 JacquetteL3-486
H3/H4 JacquetteL3-487
14/15 Posterior Jacquette

Nevi SCALERS

L3-488 Nevi 1
Anterior

An elongated disk end for anterior lingual stain and calculus removal paired with a super thin, anterior curved sickle.

L3-489 Nevi 1/H5
Anterior

An elongated disk end for anterior lingual stain and calculus removal paired with a larger H5 sickle – a thin curved blade for interproximal and cervical anterior areas.

L3-490 Nevi 2
Posterior

A super thin posterior curved sickle scaler. A unique contra-angle design great for interproximal reach and ergonomic hand positioning.

L3-491 Nevi 3
PosteriorL3-492 Nevi 4
Posterior



SCALERS

LEXER
International




Scaler for manual removal
of over-gingival backside tartar



Goldman-Fox
universal curette

SICKLE SCALERS



L3-495
U15
Towner

L3-493
H6/H7
Hygienist

L3-494
6/7
Goldman

L3-496
H5/U15
Hygienist-Towner

L3-497
U15/30
Towner-Jacquette

L3-498
H5/33
Hygienist-Jacquette

L3-499
U15/33
Towner-Jacquette


PREMIUM QUALITY DENTAL INSTRUMENTS

CURETTES



- Кюрети
- Curetas
- Küretten
- كرتسن

Alveolar bone curette for after-extraction cleaning



Hemingway alveolar bone curette for after-extraction cleaning

Hemingway



L3-500

Fig.0



L3-501

Fig.4



L3-502

Fig.5



L3-503

Fig.6

Lucas



L3-504

Fig.84 Lucas



L3-505

Fig.85 Lucas



L3-506

Fig.86 Lucas



L3-507

Fig.87 Lucas



L3-508

Fig.88 Lucas




P.K. THOMAS WAXING INSTRUMENTS

P.K. THOMAS WAXING INSTRUMENTS

Aluminium Handle
P.K. Thomas code color light handle 8mm
P.K. THOMAS WAXING INSTRUMENTS


PREMIUM QUALITY DENTAL INSTRUMENTS

L3-509
PKT-1L3-510
PKT-2L3-511
PKT-3L3-512
PKT-4L3-513
PKT-5L3-514
PKT-1L3-515
PKT-2L3-516
PKT-3L3-517
PKT-4L3-518
PKT-5



- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

P.K. THOMAS WAXING INSTRUMENTS

L3-519
PKT-1L3-520
PKT-2L3-521
PKT-3L3-522
PKT-4L3-523
PKT-5



CEMENT SPATULAS

- Шпатель для цемента
- Spatules à ciment
- Zementspatel
- Шпатель для цементу
- Espátulas para cemento
- ملاعق الأسمنت




L3-524
Fig. 1 Cement Spatulas



L3-525
Fig. 2 Cement Spatulas




L3-526
Fig. 3 Cement Spatulas



L3-527
Fig. 00 Cement Spatulas



L3-528
Fig. 0 Cement Spatulas



L3-529
Fig. 1 Cement Spatulas



L3-530
Fig. 2 Cement Spatulas



L3-531
Fig. 3 Cement Spatulas



L3-532
Fig. 4 Cement Spatulas



L3-533
Weston Cement Spatulas




L3-534
Cement Spatulas



L3-535
Cement Spatulas



Wax and Modelling Carvers

L3-536
Le CronL3-537
ZahleL3-538
ZahleL3-539
ZahleL3-540
GritmanL3-541
BealeL3-542
RoachL3-543
Vehe



ARTICULATORS

- Артикуляторы
- Articular
- Artikulations
- Артикулятори
- Articuler
- أرتيكولاتورس



L3-544

Technical Data:

Intercondylar distance	120 mm.
condylar inclination	fixed at 30°
Bennet angle	fixed at 7°
Incisal table inclination	10°
Distance between mounting plates	85 mm.
Weight	400 gr.
Opening (item 5000)	125°
Opening (item 5010)	180°



L3-545

Technical Data:

Intercondylar distance	110 mm.
condylar inclination	fixed at 30°
Bennet angle	fixed at 7°
Incisal table inclination	flat
Distance between mounting plates	85 mm.
Weight	450 gr.
Opening	180°



L3-546

Stainless Steel Finishing

L3-547

Royal Blue



L3-548



L3-549



L3-550



L3-551



L3-552



L3-553



L3-554



L3-555



L3-556



L3-557
Nakawalni
Dental Laboratory
Instrument



ALGINATE, PLASTER SPATULAS, WAX KNIVES PLASTER FOR KNIVES & PLASTER SHEAR

PREMIUM QUALITY DENTAL INSTRUMENTS



L3-558
Wooden Handle
Fig. 5
21.5 cm



L3-559
Wooden Handle
Fig. 7
19.5 cm curved



L3-560
Wooden Handle
Fig. 8
20 cm curved



L3-561
Wooden Handle
Fig. 9
22.5 cm curved



L3-562
Hook Knives



L3-563
Fig. 1
Lessmann
12.5 cm



L3-564
Fig. 2
Lessmann
17 cm



L3-565
Fig. 1
with small spoon
12.5 cm



L3-566
Fig. 2
with large spoon
17 cm



L3-567
Bacolite Handle
Wax Knives,
Small



L3-568
Bacolite Handle
Wax Knives,
Large



L3-569
Wooden Handle
Wax Knives,
Small



L3-570
Wooden Handle
Wax Knives,
Large



L3-571
Alginate
Spatulas



L3-572
Plaster Knife



L3-573
Fig. 1
17 cm



L3-574
Fig. 3
16 cm



L3-575
Fig. 4
16 cm



L3-576
Special Inox
19.5 cm



L3-577
Metal tip
220 mm
Hammer With Wooden Handle



L3-578
Metal tip
220 mm



L3-579
1
For plaster



L3-580
2
For metal



L3-581
Fret saw
250 mm



L3-582
Small
Saw Blades




L3-583
Large
Saw Blades



L3-584
Ceri Saw



L3-585
Plaster Shears
200 mm



L3-586 250 mm
L3-587 320 mm



L3-588
Ceri Saw
180 mm



L3-589
Blade
0,12x5 mm
Conf. 12 pcs



L3-590
Blade
0,25x5 mm
Conf. 12 pcs



L3-591
Blade
0,25x3 mm
Conf. 12 pcs

L3-592 Adult Palatal /
L3-593 Child Palatal



L3-594 Buccal

L3-595 Adult Palatal /
L3-596 Child PalatalL3-597 Adult Palatal /
L3-598 Adult BuccalL3-599 Lingual upper /
L3-600 Lingual LowerL3-601 Adult Buccal /
L3-602 Child BuccalL3-603 Child Palatal /
L3-604 Child BuccalL3-605
Palatal-AdultL3-606
Palatal-Adult LargeL3-607
Palatal-ChildL3-608
Buccal-SmallL3-609
BuccalL3-610
LingualL3-612
Mirror with Electric
Battery



MAXILLOFACIAL SURGERY INSTRUMENTS

LEXER
International




PREMIUM QUALITY DENTAL INSTRUMENTS



Maxiofacial Surgery Instruments

A Bard-Parker handle with a No. 15 blade is the usual scalpel used to make incisions on the gums where a flap can be raised to gain access to tooth.



Retractors

Tongue, cheek and flap retractors are needed to provide adequate access to the tooth and are legion in number.

Needle Holders and Sutures

Sutures are used to replace and secure the gum flap in its former position, to limit the size of the wound and therefore restrict subsequent contamination with food debris, and to assist in controlling bleeding. The usual suture size for intra-oral work is 4/0 and the material may be absorbable or non-absorbable Scissors Dissecting scissors and suture-cutting scissors are usually used.



Artery Forceps

These are designed to hold soft tissue without damaging it. College tweezers are used to lift up sutures prior to removing them.




Periosteal Elevators

A number are available and are similarly used like a normal elevator but periosteal elevators have two heads – one at each end. One end is sharp-tipped to expand the tooth socket while the other has a blunt end to remove granulation tissues or raise a gum flap.

Chisels and Mallet

Chisels are more appropriate in the operating theatre and are particularly useful for disto-angular third molars and upper wisdom teeth extraction.



Bone Curettes

The Mitchells trimmer is probably the most valuable instrument in the category and is used to remove granulation tissue in the tooth socket.


Bone Ronguer and Bone File

Used to smoothen any sharp bone edges after tooth extraction is done.



SCALPEL HANDLES

- Ручки для скальпелей
- Manches de bistouris
- Skalpellgriffe
- Ручки для скальпелей
- Mangos de bisturi
- مشرط مقابض



Stainless steel scalpel handle

Scalpel Handles



Locks At 90°, 60°, 30°,
Straight, 120°, 150°, 180°,
Angles



L4-001
Fig. 3
125mm

L4-002
Fig. 3
125mm

L4-003
Fig. 7
165mm

L4-004
Fig. 5
145mm

L4-005
Fig. 5A
145mm

L4-006
6 Way Adjustable
Angled
Scalpel Handle

L4-007
Handle for
Micro-Blade
13.5cm

L4-008
Fig. 4
135mm

L4-009
Double Handle
No.4 + No.3
163mm

RETRACTORS - MOUTH GAGS



- Роторасширители
- Ecarteurs - Ouvre-bouches
- Vestibulumabhalter / Mundspanner
- Роторозширювачі
- Separadores - Abre-bocas
- الكامشات-الكمامات الفم



COLUMBIA
L4-010
14cm



STERNBERG
L4-011
16cm



L4-012
Oringer
Selfretaining
cheek-lip retractor
small

L4-013
Oringer
Medium

L4-014
Oringer
Large



L4-015
Wangen Holder
150mm, with or without
lip protection



L4-016
Wangen Holder
150mm, with
or without lip protection



L4-017
Wangen Holder
Wangen holder with dual
function, 170mm



RETRACTORS - MOUTH GAGS

- Роторасширители
- Ecarteurs - Ouvre-bouches
- Vestibulumabhalter / Mundspanner
- Роторозширювачі
- Separadores - Abre-bocas
- الكامشات-الكمامات الفم



Mouth Gag Molt
L4-018 90mm
L4-019 100mm
L4-020 110mm
for (Children)

L4-021
Mouth Gag Molt
140mm for (Adults)

L4-022 105mm
L4-023 135mm
Weitlaner
Fig. 1

L4-024 105mm
L4-025 135mm
Weitlaner
Fig. 2



L4-026
Stainless Steel
Double
Ended Cheek
Retractors

L4-027
Stainless Steel
Single
Ended Cheek
Retractors

L4-028
Denhard
120mm

L4-029
Roser-Konig
160mm

L4-030
Roser-Konig
190mm

RETRACTORS



- Ранорасширители
- Ecarteurs
- Wundhaken
- Ранорозширювачі
- Separadores
- الكامشات

CHEEK RETRACTORS / TISSUE RETRACTORS

Used to hold mucoperiosteal flaps, cheeks, lips, and tongue away from the surgical area.

AUSTIN CHEEK RETRACTORS



L4-031
Fork



L4-033
Halfmoon



L4-032
Regular



L4-034
Tissue
retractor
with
handle



L4-035
Farabeuf Retractor
12cm



L4-036
Hartmann Retractor
15cm 72x24mm



L4-037
University of minnesota
Retractor 14cm



L4-038
Small 14cm

L4-039
Large 15cm

Weider cheek &
tongue depressor



RETRACTORS

- Ранорасширители
- Ecarteurs
- Wundhaken
- Ранорозширювачі
- Separadores
- الكامشات



L4-040



56X11



L4-042



55X11



L4-041
Kocher
21.5cm

L4-043
Kocher
21cm



L4-044
Cheek Retractor
Retractor and Tongue
Pushers Bishop with
Comfortable Margin,
160mm



L4-045
Cheek Retractor
Retractor Combination
250mm



L4-046
Brunings
Cheek Retractor
190mm



L4-047
Mini Retractor
About 80mm

RETRACTORS



- Ранорасширители
- Ecarteurs
- Wundhaken
- Ранорозширювачі
- Separadores
- الكامشات

15X5

8X4

L4-048
Ragnell - DavisL4-049
Senn-miller
Sharp
16cmL4-050
Senn-miller
Blunt
16cmL4-051
Appuntito
Sharp
16.5cmL4-052
Appuntito
Blunt
16.5cm

30X10

40X10

40X13

L4-056
Langenbeck
21.5cm
50X10

14X17

20X22

26X30

Middeldorf
19cmL4-060
Ramus
21.5cm
35X11



SCISSORS



L4-061
Iris
115mm, Straight



L4-062
Iris
115mm, Curved



L4-063
Micro Iris
90mm, Curved



L4-064
Iris
115mm
Angled



Goldman-Fox gum scissors
with tungsten carbide
inserts, straight



L4-065
Iris
115mm
Curved



L4-066
Goldman-Fox
130mm
Saw edge,
straight



L4-067
Goldman-Fox
130mm
Saw edge,
curved



L4-068
Goldman-Fox
130mm
Saw edge,
double curved




L4-069
La Grange
115mm
Saw edge,
double curved



L4-070
Neumann
120mm
Saw edge,
double curved



L4-071
Kilner (Ragnell)
120mm
Straight



L4-072
Kilner (Ragnell)
120mm
Curved



L4-073
Joseph
140mm
Straight



L4-074
Joseph
140mm
Curved



L4-075
Kelly
160mm
Straight




L4-076
Kelly
160mm
Curved



SCISSORS



Spencer ligature scissors
for deep suture



L4-077
Locklin
160mm
Saw edge,
angular

L4-078 120mm
L4-079 160mm
Saw edge,
angular

L4-080
Dean
170mm
Saw edge,
angular

L4-081
Spencer
90mm

L4-082
Spencer
130mm

L4-083
Spencer
115mm



L4-084
130mm
Straight

Mayo
L4-085 130mm L4-088 130mm
L4-086 160mm L4-089 160mm
L4-087 180mm L4-090 180mm
Straight Curved

Mayo
L4-091 140mm L4-093 140mm
L4-092 200mm L4-094 200mm
Straight Curved

Metzenbaum
L4-095 145mm L4-097 145mm
L4-096 200mm L4-098 200mm
Straight Curved



L4-099
Iris
115mm
Straight



L4-100
Iris
115mm
Curved



L4-101
Goldman-Fox
130mm
Saw edge, straight




L4-102
Goldman-Fox
130mm
Saw edge, Curved



Mayo
L4-103 140mm
L4-104 170mm
L4-105 230mm
Curved




Metzenbaum
L4-106 145mm
L4-107 180mm
L4-108 200mm
Straight



Metzenbaum
L4-109 145mm
L4-110 180mm
L4-111 200mm
Curved



L4-112
Kelly
160mm
Straight



L4-113
Kelly
160mm
Curved




SCISSORS





NEEDLE HOLDERS





SCISSORS

LEXER
International



Noyes/Castroviejo
gum scissors
with curved and
very thin beaks



L4-137
Castroviejo
Straight 12cm



L4-138
Castroviejo
Curved 12cm



L4-139
Castroviejo
Straight 14cm



L4-140
Castroviejo
Curved 14cm



L4-141
Micro Castroviejo
Curved 18cm



L4-142
Westcott
Curved 11cm



L4-143
Noyes
110mm, Straight



L4-144
Noyes
110mm
Curved



L4-145
Noyes
110mm
Angular



L4-146
Castroviejo
150mm
Curved

NEEDLE HOLDERS



- Иглодержатели
- Porte-aiguilles
- Nadelhalter
- Голкотримачі
- Porta-agujas
- أصحاب إبرة



L4-147
Delicate
150mm



Adson
L4-148
170mm



Mayo-Hegar
L4-149 130mm
L4-150 140mm
L4-151 160mm
L4-152 180mm
L4-153 200mm



L4-154
Crile-Wood
150mm



Olsen-Hegar Needle
Holder With Scissors
L4-155 140mm



Olsen-Hegar needle holder
and scissors combined



Gillies Needle
Holder With Scissors
L4-156 160mm



NEEDLE HOLDERS

- Иглодержатели
- Porte-aiguilles
- Nadelhalter

- Голкотримачі
- Porta-agujas
- أصحاب إبرة



L4-157
Delicate
150mm




Mayo-Hegar
L4-158 130mm
L4-159 140mm
L4-160 160mm
L4-161 180mm
L4-162 200mm



L4-163
Crile-Wood
150mm
straight



L4-164
Crile-Wood
150mm
curved



Mini-Ryder needle holder with tungsten carbide inserts



Crile
L4-165 130mm
L4-166 150mm
L4-167 180mm



L4-168
Halsey
130mm



L4-169
Baumgartner
130mm



Derk
L4-170 120mm
L4-171 140mm

NEEDLE HOLDERS



- Иглодержатели
- Porte-aiguilles
- Nadelhalter
- Голкотримачі
- Porta-agujas
- أصحاب إبرة



Mathieu
L4-172 140mm
L4-173 170mm



Castroviejo
L4-174 140mm Str.
L4-175 180mm Str.
L4-176 140mm Cvd.
L4-177 180mm Cvd.



Olsen-Hegar Needle
Holder With Scissors
L4-178 140mm



Gillies Needle
Holder With Scissors
L4-179 160mm



Castroviejo
L4-180 140mm Str.
L4-181 180mm Str.
L4-182 140mm Cvd.
L4-183 180mm Cvd.



L4-184
Zarral
with lock
200 mm



L4-185
Castroviejo TC
14cm



L4-186
Castroviejo TC
14cm




L4-187
Castroviejo TC
14cm



HAEMOSTATIC FORCEPS & TOWEL CLAMPS

- Зажимы артериальные, зажимы для салфеток
- Pinces hemostatiques, Pinces a compresses
- Arterienklemmen, Tuychklemmen
- Затискачі артеріальні, затискачі для серветок
- Pinzas hemostáticas, Pinces para compresas
- الملفت هايموستاتيك المشابك منشفة





L4-224
curved



Foerster
Dressing Forceps
L4-225 200mm



Young
Tongue Forceps
L4-226 170mm



Collins
Tongue Forceps
L4-227 170mm



Spare rubber



L4-228 20 cm
L4-229 25 cm
L4-230 30 cm



Rogge
L4-231
27 cm





Cheate
Box Joint
L4-232 27 cm



RUBBER DAM INSTRUMENTS

LEXER
International



PREMIUM QUALITY DENTAL INSTRUMENTS



- Резиновая дамба зажимы комплекты
- Trousses de caoutchouc
- Gummidammklammern
- Гумова дамба затискачі комплекти
- Conjuntos de abrazaderas de retención de caucho
- المطاط السد المتشابك مجموعات

Rubber Dam Kits

**L4-246**

On request the Content of the rubber dam kit can be changed with other instruments.

Content:

- 01 Anisworth
- 02 Brewer
- 03 Clamps Stand
- 04 Rubber Dam Frame L
- 05 Rubber Dam Frame M
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case

08 Rubber Dam Clamp Fig. 212

09 Rubber Dam Clamp Fig. 211

10 Rubber Dam Clamp Fig. 8

11 Rubber Dam Clamp Fig. 8A

12 Rubber Dam Clamp Fig. 14

13 Rubber Dam Clamp Fig. 14A


14 Rubber Dam Clamp Fig. 0

15 Rubber Dam Clamp Fig. 00

16 Rubber Dam Clamp Fig. 2

17 Rubber Dam Clamp Fig. 2A

18 Rubber Dam Clamp Fig. 7

**L4-247**

On request the Content of the rubber dam kit can be changed with other instruments.

Content:

- 01 Ivory
- 02 Ivory light weight
- 03 Clamps Stand
- 04 Rubber Dam Frame M
- 05 Rubber Dam Frame L
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case

08 Rubber Dam Clamp Fig. 209

09 Rubber Dam Clamp Fig. 212

10 Rubber Dam Clamp Fig. 211

11 Rubber Dam Clamp Fig. 8

12 Rubber Dam Clamp Fig. 8A

13 Rubber Dam Clamp Fig. 14

14 Rubber Dam Clamp Fig. 14A

15 Rubber Dam Clamp Fig. 0

16 Rubber Dam Clamp Fig. 00

17 Rubber Dam Clamp Fig. 2

18 Rubber Dam Clamp Fig. 2A

19 Rubber Dam Clamp Fig. 7



RUBBER DAM CLAMPS

- Кламмера для коффердама
- Crampons pour digue
- Kofferdam Klammern
- Кламмера для коффердама
- Grapas para dique de caucho
- المشابك المسد المطاط



L4-248
Fig. 15



L4-249
Fig. 18



L4-250
Fig. 24



L4-251
Fig. 25



L4-252
Fig. 26



L4-253
Fig. 28



L4-254
Fig. 30



L4-255
Fig. 31



L4-256
Fig. 51



L4-257
Fig. 138



L4-258
Fig. 139



L4-259
Fig. 200



L4-260
Fig. 201



L4-261
Fig. 202



L4-262
Fig. 203



L4-263
Fig. 204



L4-264
Fig. 205



L4-265
Fig. 3



L4-266
Fig. W3



L4-267
Fig. 4



L4-268
Fig. 56A



L4-269
Fig. 7



L4-270
Fig. 7A



L4-271
Fig. 8



L4-272
Fig. 8A



L4-273
Fig. W8A



L4-274
Fig. 10



L4-275
Fig. 11



L4-276
Fig. 12A



L4-277
Fig. 13A



L4-278
Fig. 14



L4-279
Fig. 14A



L4-280
Fig. 14AN



L4-281
Fig. 22



L4-282
Fig. 27



L4-283
Fig. 207



L4-284
Fig. 29



L4-285
Fig. 206



L4-286
Fig. 208



L4-287
Fig. 209



L4-288
Fig. 1



L4-289
Fig. 1A



L4-290
Fig. 2



L4-291
Fig. 2A



L4-292
Fig. W2A



L4-293
Fig. 0



L4-294
Fig. 00



L4-295
Fig. 210



L4-296
Fig. 211



L4-297
Fig. 212



L4-298
Fig. 6



L4-299
Fig. 9



L4-300
Fig. N06



L4-301
Fig. U67



L4-302
Fig. P1



L4-303
Fig. P2



L4-304
Fig. G-1



L4-305
Fig. G-2

CROWN REMOVER



- Инструменты для удаления коронок
- Arrache-couronnes
- Teleskopkronenzangen
- Інструменти для видалення коронок
- Arranca-corona
- مزيلات التاج




L4-306
Crown Spreaders




L4-307
Mershon
140mm
Band Pusher




L4-308
Crown Spreader



L4-309
Crown Remover



L4-310
Crown Remover



L4-311
Crown Remover



CROWN REMOVER

- Инструменты для удаления коронок
- Arrache-couronnes
- Teleskopkronenzangen
- Інструменти для видалення коронок
- Arranca-corona
- مزيلات التاج



L4-312
Turkey Pattern



L4-313
Manual Pattern




L4-314
Crown Opener



L4-315
Temporary
Crown Removing Plier



- Инструменты для удаления коронок
- Arrache-couronnes
- Teleskopkronenzangen
- Інструменти для видалення коронок
- Arranca-corona
- مزيلات التاج



Terms of Use

Removal of:-

- Permanently cemented crowns and /or bridges
- Temporarily cemented provisional or inal restorations
- Maryland bridges
- Orthodontic brackets/bands
- Removable appliances "stuck" at attachments
- Multi-unit restoration "stuck" on working lab model

Basic Kit Includes:-

- Crown and Bridge Remover # L4-316 / L4-317
- Single crown remover tip # A
- Crown remover tip Biid # B
- Crown & bridge removing tip Sickle # C
- Single crown remover Contra angle tip # D
- Single crown remover tip Inverted curved # E
- Single Crown remover tip straight # F

- Bridge remover Flexi-wire set # L4-318
- Bridge remover tip with lock # L4-319
- Tips adjustment Key # G
- Open-end adjustment wrench # H

Precautions

1. To remove crowns placed over pins or screw-type posts, use lowest "strike intensity" to prevent fracture.
2. This instrument could cause damage to the margins of soft gold crowns. If re-use is anticipated, cut a small groove parallel to the lingual margin.
3. Highest strike intensity could cause root fracture; be careful in long span restorations.
4. Only low intensity strike force should be used with teeth having high mobility and little bone support.
5. For highest strike intensity settings, let the individual situation and common sense prevail.



MATRIX RETAINERS

LEXER
International

- Матрицедержатели
- Porte Matrices
- MatrizenSpanner
- Матрицеутримувачи
- Porta Matrices
- مصافحة الخدم



L4-320
Fig. 5
Ivory 8
5 cm



L4-321
Fig. 6
Ivory 8
6 cm



L4-322
Fig. 7
Ivory 8
7 cm



L4-323
Fig. 1
Tofflemire-Universal

Siqveland Retainers



L4-324
Fig. 2
Tofflemire-Junior



L4-325
Fig. 3
Tofflemire-Senior



L4-326
Narrow



L4-327
Wide



L4-328
Stainless Steel
Amalgam Well



L4-329
Ivory 12



L4-330
Ivory 1



L4-331
Ivory



L4-332
Elliot

AMALGAM INSTRUMENTS





- Инструменты для амальгамы
- Instruments pour amalgame
- Amalgam instruments
- Інструменти для амальгами
- Instrumentos para amalgama
- الصكوك الملغمة

AMALGAM CARRIERS



TIP SIZES FOR DOUBLE-END





AMALGAM INSTRUMENTS

LEXER
International

- Инструменты для амальгамы
- Instruments pour amalgame
- Instrumentos para amalgama
- Amalgam instruments
- المكواكب الملحمة





- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente fur die parodontologie
- Інструменти пародонтологічні
- Instrumentos para periodoncia
- الصكوك اللثة

PERIOSTEALS ELEVATORS

Molt universal periosteal elevator M8

L4-381
Fig. M7
Molt

L4-382
Fig. M8
Molt

L4-383
Fig. M9
Molt

L4-384
Fig. 2
Hourgen

L4-385
Fig. PR3
Prichard

L4-386
Fig. 9A
Allen


L4-387
Fig. GF
Goldman-Fox

L4-388
Williger
16cm



BONE FILE

LEXER
International



PREMIUM QUALITY DENTAL INSTRUMENTS



Miller



L4-406
Fig. 22
Seldin
19cm



L4-407
Fig. 23
Seldin
19cm



L4-408
Fig. 2
West



L4-409
Fig. 9
Molt



L4-410
Fig. 14
Goldman Fox



L4-411
Prichard



L4-412
Fig. 1
Freer



L4-413
Fig.2
Freer



L4-414
Fig. 1
Williger



L4-417
Fig. 4



PERIOSTEAL ELEVATORS - MALLETS & BONE CHISEL

MALLETS

Used with a chisel to split teeth and/
remove bone



L4-418
Mead1
19cm



L4-419
Mead2
17cm

BONE CHISEL

Used to scoop bone and make trough



L4-420
Gardner
3mm



L4-421
Gardner
4mm



L4-422
Chandler
5mm



L4-423
Chandler
7mm



L4-424
Partsch
19cm

L4-425
Mead
19cm



L4-426
Plaster inserts
Autoclaveable max



L4-427
Fig. 1
Freer



L4-428
Fig. 2
Freer



L4-429
Fig. 3
Freer



L4-430
Fig. 1
3.8mm
Straight



L4-431
Fig. 2
3.8mm
Curved



L4-432
Fig. 3
7.5mm
Curved




L4-433
Fig. 4
7.5mm
Curved

PERIODONTAL INSTRUMENTS



- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente fur die parodontologie
- Інструменти пародонтологічні
- Instrumentos para periodoncia
- الصكوك للثنا



L4-434
Wakefield 1W
170mm

L4-435
Wakefield 2W
170mm

Lucas 16cm
L4-436 4mm
L4-437 6mm

Lucas 16cm
L4-438 3mm
L4-439 5mm

Partsch 140 mm
L4-440 3mm
L4-441 4mm
L4-442 5mm
L4-443 6mm

Partsch 170 mm
L4-444 3mm
L4-445 4mm
L4-446 5mm
L4-447 6mm

Partsch 140 mm
L4-448 3mm
L4-449 4mm
L4-450 5mm
L4-451 6mm

Partsch 170 mm
L4-452 3mm
L4-453 4mm
L4-454 5mm
L4-455 6mm

Osteotome Chisels - Set of 4



L4-456
1
Ochsenbin

L4-457
2
Ochsenbin

L4-458
7.5 mm Straight
Osteotome
(7-10-13-15-18 mm markings)

L4-459
4 mm Straight
Osteotome
(7-10-13-15-18 mm markings)


L4-460
7.5 mm Curved
Osteotome
(7-10-13-15-18 mm markings)

L4-461
4 mm Curved
Osteotome
(7-10-13-15-18 mm markings)



PERIODONTAL INSTRUMENTS

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente fur die parodontologie
- Інструменти пародонтологічні
- Instrumentos para periodoncia
- الصكوك اللثة



BONE ROOUNGER

- Костные кусачки
- Hohlmeißelzangen
- Pinces-gouges
- Кусачки для кісток
- Pinzas-gubias
- رونجیورس العظام

BONE RONGEURS

Used for trimming and recontouring alveolar bone and gross tissue removal



Blumenthal L4-477



Blumenthal 45° L4-478



Blumenthal 90° L4-479

L4-483
Double Action

Mead L4-480



Friedman L4-481




Mini Friedman L4-482

L4-484
BlumenthalL4-485
Blumenthal 45°L4-486
Blumenthal 90°L4-487
Micro Friedman
with La Grange Bending
(s-shape)




FILM HANGER, NAPKIN HOLDER & MOUTH PROBS

LEXER
International



Key Rings







ORTHODONTIC PLIERS & INSTRUMENTS

orthodontics pliers in Clinical Dentistry Orthodontics and Dentofacial Orthopaedics

Driven by the demand for tooth alignment, general practitioners are increasingly being asked to provide orthodontic treatment.


Dentists keen to participate in this area of expertise need to demonstrate the skills to identify suitable patients, the ability to provide conventional and contemporary orthodontic treatment and be aware of its limitations.

This is exactly what this programme delivers, along with the vital core competencies to excel in this unique postgraduate area.

The objective of the Orthodontics and Dentofacial Orthopaedics is to provide dental surgeons with the core knowledge and clinical skills essential for the delivery of competent, high-quality dental care.

It encompasses a blended learning approach and consists of a combination of lectures, seminars, workshops, clinical case discussions and laboratory sessions.

The flexible modular structure enables the programme to meet a wide range of student needs.





- Щипцы для ортодонтии и протезирования
- Pince pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- كماشة لقصويم الأسنان بروشتيتس



L5-001
Cap Removing Plier



L5-002
1-Distal End Cutter
This instrument cut wires close to the buccal tube and safely holds the distal end. The tungsten Carbide insert ensure precision cutting, long wear and unsurpassed cutting performances.
Maximum Cutting Capacity.
arch wires up to .021"x.025"(.55mx.64mm)



L5-003
1-Distal End Cutter Flush Cut with Safety Hold
This instrument cut wires close to the buccal tube and safely holds the distal end. The tungsten Carbide insert ensure precision cutting, long wear and unsurpassed cutting performances.
Maximum Cutting Capacity.
arch wires up to .021"x.025"(.55mx.64mm)



L5-004
Micro Ligature Cutter
Maximum Cutting Capacity.
.008" (20mm)
Maximum Cutting Capacity.
.012" (30mm)



L5-005
1-Distal End Cutter Flush Cut
This instrument cut wires close to the buccal tube and safely holds the distal end. The tungsten Carbide insert ensure precision cutting, long wear and unsurpassed cutting performances.
Maximum Cutting Capacity.
arch wires up to .021"x.025"(.55mx.64mm)



L5-006
Hard wire cutter 15° Degree Angle
Easily Cuts any Size arch wire or Pin.
The inserted tips ensure precision cutting and perfect tip alignment.
Minimum Cutting and perfect tip alignment. Minimum Cutting Capacity.
.012"(30mm)
Maximum Cutting Capacity:
arch wires up to 012"x025"
(.55mmx.64mm)



L5-007
Hard wire cutter Standard Straight
Easily Cuts any Size arch wire or Pin.
The inserted tips ensure precision cutting and perfect tip alignment.
Minimum Cutting and perfect tip alignment. Minimum Cutting Capacity.
.012"(30mm)
Maximum Cutting Capacity:
arch wires up to 012"x025"
(.55mmx.64mm)



L5-008
Hard wire cutter, easily cuts any size arch wire or pin. The inserted tips ensure precision cutting and perfect tip alignment.
Minimum Cutting Capacity: .012".(30mm)
Maximum Cutting Capacity:
arch wires up to .021"x.025"(.55mmx.64mm)



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pince pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- كماشة لقصويم الأسنان بروسيتيس



L5-009
Pin Ligature cutter 15° degree angled tips to allow easy access and better visibility to reach in Difficult areas. Inserted tips ensure precision cutting of Ligature wire, pins and elastics.
Minimum Cutting Capacity: .008".(20mm)
Minimum Cutting Capacity: .016".(41mm)



L5-010
Standard Ligature cutter
This cutter has fine tips for easier access to reach difficult areas cuts pins, wire ligatures and elastics.
Minimum Cutting Capacity: .008".(20mm)
Minimum Cutting Capacity: .016".(41mm)



L5-011
Mini Ligature cutter 50° degree angled upper the narrow tip is angled at 50° degree for excellent access to the upper arch. The inserted tip ensure precision cutting of ligature wires and elastics.
Minimum Cutting Capacity: .008".(20mm)
Minimum Cutting Capacity: .012".(30mm)



L5-012
Bracket Removing Plier Angled TC - inserts
Sharp beak wedge can easily go under bracket base. Lifts off bracket from tooth surface with no discomfort to patient. Firm grip assures, easy in horizontal and vertical use.



L5-013
Ligature wire forming plier
for soft wires up to .014" (.36mm)



L5-014
Crimpable Hook Plier TC inserts
The head of these pliers is designed with a slot to fit over the base of the hook for a secure hold. One squeeze and the hook is secured on the arch wire.



L5-015
Cinch back plier
The Serrated tips of this instrument hold the wire securely while bending the wire 90° degree distal to the buccal tube, awoinding patient for use with Nickle Titanium wires .012"x.018"(.30mmx46mm



L5-016
Johnson contouring plier
Design to add contour or Curvature by crimping trail-fitted Crowns. This Contouring / crimping ensures that the crows fits securely on the patients tooth. slim weingart pliers



L5-017
Adams Plier 325 with hole for retaining elements



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- كماشة لتقويم الأسنان بروشتيتيس

**L5-018****How Plier (straight)**

A versatile utility plier with a long tapered tips and have serrated tips (Also available) in TC-Tip to prevent wire slippage.

Point Size: 2mm,4mm
Also Available in TC-Tip

**L5-019****How Plier (Cvd)**

Similar to How Plier (straight) but with 45° degree angled tips for convenient access to lingual and to reach difficult areas. Matching serrated tips hold wires firmly at a convenient working angle.

Point Size: 2mm,4mm
Also Available in TC-Tip

**L5-020****Direct Bond Remover (Angled)**

Some Features as Direct Bond Remover plier (straight) except the 60° degree angle allows easier access in those to reach difficult posterior areas.

**L5-021****Direct Bond Remover (straight)**

Force is applied to the bracket tooth interface. The tips are positioned mesially, distally or at the occlusal gingival interface to remove direct bond brackets Safely and Comfortably.

**L5-022****Posterior Band Remover (Long)**

Some features as posterior band remover plier (short) but with a longer lower band removing arm. Replacement pads are available in 3/16" and 1/8" (4mm) diameters.

Also Available in TC-Tip

**L5-023****Posterior Band Remover (Short)**

The plastic pad rests on the occlusal surface of the molar which positions to the lower removing tip in the perfect position under the gingival edge of the band. Replacement pads are available in 3/16" and 1/8" (4mm) diameters.

Also Available in TC-Tip

**L5-024****Weingard Plier Slim**

Precision serrated tips are designed for heavy wires and to prevent the wire from slipping. For easy and safe insertion or removed of all arches.

Holds wire up to .030"(.76mm)
wire bending capacity: .030"(.76mm)

Also Available in TC-Tip

**L5-025****Weingard Plier Utility**

Carefully tapered beaks fit easily between brackets and are rounded for safety and comfort. Holds round, Square and rectangular wires up to .22"x.028"(.56x.71mm)

Also Available in TC-Tip

**L5-026****Torque Bending Plier Pair use**

Designed to place a torque in a specific section of arch wire with out distortion in the reminder of the wire. The pliers use a key (either .016"/.018" or .018"/.022") that torque the section between holding beaks. Labial and lingual torque can be Placed Adjacent to each other easily. Torque can be limited to single tooth.



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pince pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- كماشة لتقويم الأسنان بروسثيتس



L5-027

Tweed plier slim

Designed for placing first, second and third order bends wires up to .022"x.025"(.53mmx.64mm)

Also Available in TC-Tip



L5-028

Lingual Arch Bending plier

Grooved to hold Several Sizes of wire, Jaw Thickness is .12"(3.05mm) I deal for double and triple back bends. Bends round, Square and rectangular wire up to .036"(0.9mm)max.



L5-029

Band Crimping Plier

This tips allow for easy access into molar and premolar bands.



L5-030

Tweed Arch Forming Plier

Designed for torquing and bending square or rectangular wires. The blades are designed to be parallel at .020" opening. Can be used on wires up to .021"x.025"(.55mmx.64mm)



L5-031

Jaraback Pliers

Long, tapered design with fine detail work on serrated tips allows for intricate bends and utility work. Tip sizes are .028" for round and half-moon. Wire bending capacity: .028"



L5-032

Tweed Loop Forming Pliers

The perfect plier for making accurate, omega and closing loops. can be used with round & square wire up to .021"x.025" (.55mmx.64mm)



L5-033

Young Plier


For labial, lingual, arches and Canine loops. For hard wire up to .028"(0.7mm) max. wire Bending Capacity (.028"(0.7mm) C-1280 Aderer Small (3 jaw plier Small) Same features as the 3 jaw plier but with Slightly Smaller tips for working with wire up to .020" (.51mm) wire Bending Capacity (.020"(.51mm) C-1281Aderer Precision aligned tip are rounded for consistent contouring and bending without damaging the archwire. One-piece construction fo the double tip ensures Superior strength. Can be used on wire up to .030"(.76mm) wire Bending Capacity .030"(.76mm)



L5-034

Nance Plier

Precision loop forming with four step tips 3mm,4mm,5mm and 6mm, Excellent for round and square wires up to .022" (.56mm) All working edges have been beveled to prevent wire Scoring. wire bending capacity. 022"(.56mm)



L5-035

These Multipurpose Specialty Pliers

are used for easily forming stops, tie backs, or Adams clasps for retainer.

PLIERS FOR ORTHODONTICS & PROSTHETICS



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- كماشة لتقويم الأسنان بروشتيتيس



L5-036
V Bend Plier

Creates a 1 mm v bend to any part of the archwire. Provides positive stops and prevent archwires from "walking"



L5-037
De la Rosa with three Grooves
Arch contouring and forming plier
with three grooves, for wires .016", .018" and .022"



L5-038
De la Rosa without Grooves
perfect for Subtle Contouring of arch wires up to .030" (.76mm)
wire Bending Capacity: .030" (.76mm)



L5-039
Adams plier with Hole
Sling Claps and wire bending pliers
Also Available in TC-Tip



L5-040
Slim Bird Beak Plier
This popular loop forming and bending plier has inserted tips that provide excellent wear. The precise round and square tips bend arch wires from .012" to .020" (.30mm-.51mm)
Wire bending capacity: .020" (.51mm)
Also Available in TC-Tip



L5-041
Bird Beak Pliers
The most popular and versatile plier for working with round wire up to .028" (.71mm).
The round beak is precision manufacture to a perfect cone wire bending capacity: .028" (.71mm)
Also Available in TC-Tip



L5-042
Light wire Bird Beak Plier
A longer cone than a standard bird beak is useful for forming precise loops in light wires up to .020" (.51mm)
wire Bending Capacity: .020" (.51mm)



Aderer
L5-043 11.5 cm,
L5-044 13cm



L5-045
Reynolds contouring pliers
Also Available in TC-Tip



PLIERS FOR ORTHODONTICS & PROSTHETICS

LEXER
International

- Щипцы для ортодонтии и протезирования
- Pince pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- كماشة لتقويم الأسنان بروسثيتس



Step Plier
L5-046 0.5 mm
L5-047 0.75 mm
L5-048 1 mm



L5-049
Arrow Clasp Bending Pliers



L5-050
Schwarz Pliers



L5-051
Lingual Arch Bending Pliers

ORTHODONTICS CUTTERS TC




L5-052
Adhesive Removing Plier



L5-053
Canine
Contouring



L5-054
Ball & Socket



L5-055
Hard Wire Cutter



- Щипцы для ортодонтии и протезирования
- Pince pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- كماشة لتقويم الأسنان بروسثيتس



L5-056
Optical



L5-057
Anterior Band
Remover



L5-058
Band/Crown
Crimping



L5-059
Ligature Tying



L5-060
Separating Module



L5-061
Left
Band Pinching



L5-062
Right
Band Pinching



MEASURING INSTRUMENTS



L5-063
Clamp, Osteometer
Length about 155 mm,
ridge width measurement,
two-sided scale, measuring
range 0 - 25 mm



L5-064
Measuring Loop
straight length 90 mm,
measuring range
0 - 40 mm,



L5-065
Measuring Loop
angled, length 90 mm,
measuring range
0 - 40 mm,



L5-066
Measuring loop, Castroviejo
angled, length 175 mm,
measuring range 0 - 40 mm,
with window



L5-067
Measuring Loop
angled, length 175 mm,
measuring range
0 - 40 mm,



L5-068
Measuring Loop
straight length 175 mm,
measuring range
0 - 40 mm,



L5-069
Fig. 1
Iwanson
spring caliper For metal
and porcelain.



L5-070
Fig. 2
Iwanson
spring caliper
For wax.



Caliper / Calipers Rule



L5-071
Delux Boley
Gauge



L5-072
Boley Gauge



L5-073
Beerendonk
130mm
Mini Gauge



L5-074
Zuricher Gauge
130mm



L5-075
Boley Gauge
Stainless Steel



L5-076
Dentaurum Ruler



L5-078
Compasso
De Willis



L5-079
Band seater
auto cleavable



L5-080
Band seater
auto cleavable



L5-081
Band seater
auto cleavable



L5-082
Band Pusher

L5-083
Strip Saw



ORTHODONTIC INSTRUMENTS

LEXER
International

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Інструменти ортодонтичні
- Instruments Orthodontiques
- أدوات تقويم الأسنان



L5-084
Endo Finger
Ruler Span
Measure Scale



L5-085
Endo Finger
Ruler Span
Measure Scale



L5-086
Endo Finger
Ruler Span
Measure Scale



L5-087
Endo Finger Ruler
Span Measure Scale



L5-088
Endo Finger
Ruler Span
Measure Scale



L5-089
Endo Finger
Ruler Span
Measure Scale



L5-090
Endo Finger Ruler
Span Measure Scale



Falcon bracket positioning gauges
Four different heights for bracket positioning: 3.5,
4, 4.5, & 5 mm
Two different models for .018" or .022" slot brackets
anodized aluminum, autoclavable



Bracker placer for posterior teeth
Very helpful instrument for placing D.B. brackets from bicuspids to molars, the small head offers excellent view



L5-091
1

HEIGHT GAUGE MICROMETER



L5-092
2

Also available in stainless steel



L5-093
Moveable
Height Gauge



L5-094
Positive
Postioner



L5-095
Boon's Gauge



- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Інструменти ортодонтичні
- Instruments Orthodontiques
- أدوات تقويم الأسنان



Bracket Holder Tweezers



L5-096
Direct Bracket
Holder Tweezers

L5-097
Exect Bracket
Holder Tweezers

L5-098
Direct Bonding
Bracket Holders 135mm

L5-099
Buccal Tube Houxer

L5-100
Buccal Tube Houxer

L5-101
Exect Bracket
Holder Tweezers



L5-102
Band seating &
removing instruments

L5-103
Ortho Light double ended
ligature director scaler

L5-104
Ortho Light Band
Pusher Scaler

L5-105
Distal Bender
Silicon Grip

L5-106
Distal Bender
Silicon Grip

L5-107
Distal Bender
Silicon Grip



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- كماشة لتقويم الأسنان بروسثيتس



L5-108
Weingart TC



L5-109
Right



L5-110
Left



L5-110A
Ligature Plier Mathieu



L5-110B
Elastic Placer Mathieu

BAND FORMING



L5-111
#1
LIGATURE
TUCKER



L5-113
#1
UTILITY
INSTRUMENT



L5-112
#1
LIGATURE DIRECTIONAL
INSTRUMENTS



L5-114
#1
ELASTIC LIGATION
APPLICATOR



STEP 1 A hole is drilled for the implant. A cover screw is put in and left till the bone fuses and is ready for the implant.

STEP 2 The cover screw is removed in readiness for the next stage.

STEP 3 The abutment is attached to the implant.

STEP 4 The final stage where the crown (prosthetic tooth) is attached.





Instruments for dental implants are made of high quality stainless steel.

Our line of instruments include, but is not limited to, the following products:

Bone Grafting Instruments

- Bone Curettes (SPRATT) – designed for harvesting bone
- Bone Grinder – used in conjunction with a mallet to grind bone into appropriate particle sizes
- Bone Basin – weighted bottom design prevents tipping or spilling while mixing bone or other bone grafting materials
- Bone Compactors – feature serrated tips for compacting bone

Scalpels & Biopsy Punches

- Scalpels – offers convenience of a disposable with razor sharp cutting edge for superior cutting
- Biopsy Punches – designed at a steeper angle to ensure a finely honed, sharper edge
- Bone Spreading Chisels
- Freer Chisels – double or single cut tips used for widening ridges during implant preparations
- Graduated Chisels – feature reference markings for measurement during a bone ridge augmentation and are available in straight and curved patterns

Osteotomes

- Concave Osteotomes – break the floor of the maxillary sinus through the implant cavity and are available in curved (off set) and straight (with or without stops)
- Convex Osteotomes – raise the floor of the maxillary sinus and are available in curved (off set) and straight patterns (with or without stops)
- Progressive Osteotomes – designed to prepare the implant site in the maxillary area where bone density is 3 or 4
- Bone Spreading Osteotomes – sharp or blunt configurations for progressive circumferential expansion of the implant site

Sinus Lift Instruments

- Sinus Lift Curettes and Elevators – lightweight, wide diameter stainless steel hand instruments lift the sinus membrane for bone augmentation procedures

Gauges & Calipers


- Gauges & Calipers – various measuring devices available to determine appropriate length, width, and distance for implant procedures

Collagen Products

- Collagen wound dressing product offering includes products.

Sterilization Trays

- Sterilization Trays – designed to protect and store instruments during sterilization, transportation and storage





- Инструменты для имплантологии
- Implantologie
- Implantologie

- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



IMPLANT GAGES

**L6-001**

Multifunctional Meter
Length 90 mm,
different sizes,
8-10-11mm,
corrosion-resistant

L6-002

Multifunctional Meter
Length 90 mm,
different sizes,
5-13-15-17 mm,
corrosion-resistant

L6-003

Measuring Probe
Length 170 mm,
(6-8-10-14-16 mm),
corrosion-resistant

L6-004

Implant
Depth Gauge

L6-005

Implant
Depth Gauge
& Key


L6-006

Implant
Depth Gauge


L6-007

Implant
Depth Gauge

IMPLANT GAGES

**L6-008**

Double Ended Osteotomy Depth Gauge

**L6-009**

Double Ended Osteotomy Depth Gauge

**L6-010**

Double Ended Osteotomy Depth Gauge

**L6-011**

Depth & Gingival Gauges



IMPLANTOLOGY


- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان




L6-012
Ruler



L6-013
Ortho Ruler / Cone Socket



L6-014
Dynamometer



L6-015
Bone Thickness Gauge
40mm, 180mm



L6-016
Micro Boley
Gauge Twisted



L6-017
Micro Boley
Gauge Straight



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



Sinus Lift Curettes and Elevators – lightweight, wide diameter stainless steel hand instruments lift the sinus membrane for bone augmentation procedures



L6-018
Fig. 1

L6-019
Fig. 2

L6-020
Fig. 3

L6-021
Fig. 4

L6-022
Fig. 5

L6-023
Fig. 6

L6-024
Fig. 7

L6-025
Fig. 8




IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



Sinus Lift Curettes and Elevators – lightweight, wide diameter stainless steel hand instruments lift the sinus membrane for bone augmentation procedures

Periosteal Elevators



L6-026
Fig. 1

L6-027
Fig. 2

L6-028
Fig. 3

L6-029
Fig. 4

L6-030
Fig. 5

L6-031
Fig. 6




- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



Sinus Lift Curettes and Elevators – lightweight, wide diameter stainless steel hand instruments lift the sinus membrane for bone augmentation procedures

IMPLANT SINUS LIFT INSTRUMENTS HOLLOW HANDLE 9.5MM LENGTH 170MM



L6-032
Acute 6577

L6-033
Sinus1

L6-034
Sinus2

L6-035
Sinus3,
Obtuse 6578

L6-036
Sinus4

L6-037
Sinus5

L6-038
Sinus6

L6-039
Sinus7




IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



Sinus Lift Curettes and Elevators – lightweight, wide diameter stainless steel hand instruments lift the sinus membrane for bone augmentation procedures

SINUS LIFT INSTRUMENTS HOLLOW HANDLE 9.5MM LENGTH 170MM



L6-040
Sinuslift 1

L6-041
Sinuslift 2

L6-042
Sinuslift 3

L6-043
Sinuslift 4

L6-044
Sinuslift 5

L6-045
Sinuslift 6

L6-046
Membrane
3Positioning 1

L6-047
Condenser 1

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان




Sinus Lift Curettes and Elevators – lightweight, wide diameter stainless steel hand instruments lift the sinus membrane for bone augmentation procedures

Sinus Lift Instrument,
Flexible Sinus Lift instruments
reduce the risk of perforation of the
“Schneider's membrane”

Sinus Lift Elevator,
10 mm Straight grip
Length 175 mm, Sinus lift elevators, dull,
for gentle peeling and lifting the
Schneider membrane

Sinus lift instrument
10 mm Straight grip Length 180 mm,
Buser (B1, B2, B3) raspatories,
quality medical stainless steel,
corrosion-resistant



L6-048
Fig. 1

L6-049
Fig. 2

L6-050
Fig. 3

L6-051
Fig. 1

L6-052
Fig. 2

L6-053
Fig. 3

L6-054
B01

L6-055
B02

L6-056
B03




IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



Sinus Lift Curettes and Elevators – lightweight, wide diameter stainless steel hand instruments lift the sinus membrane for bone augmentation procedures

Sinus lift instrument 10 mm Straight grip Sinus lift elevators, blunt or sharply curved
(H01, H02, K03, K04, K78, K77)



L6-057
Sinus lift
instrument
H01

L6-058
Sinus lift
instrument
H02

L6-059
Sinus lift
instrument
K03

L6-060
Sinus lift
instrument
K04


L6-061
Sinus lift
instrument
K78

L6-062
Sinus lift
instrument
K77

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



Sinus Lift Curettes and Elevators – lightweight, wide diameter stainless steel hand instruments lift the sinus membrane for bone augmentation procedures





IMPLANTOLOGY


- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



Sinus Lift Curettes and Elevators – lightweight, wide diameter stainless steel hand instruments lift the sinus membrane for bone augmentation procedures

SINUS LIFT INSTRUMENT

Sinus Lift Elevators, 10 mm Hollow Handle length 170mm Blunt Curved Set of 6 Pcs



L6-070

L6-071

L6-072

L6-073

L6-074

L6-075



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Імплантологія
- زراعة الأسنان



Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm
IMPLANT SURGICAL CHISELS



L6-076 13K/TG



L6-077 13K/TG Titanium Plasma Coated Point



L6-078 13K/TG



L6-079 13KL/R Titanium Plasma Coated Point



L6-080 13KL 4/6mm



L6-081 Pt2 Titanium Plasma Coated Point
A/V Regular



L6-082 Pt1 Titanium Plasma Coated Point
A/V Regular



L6-083 Pt3 Titanium Plasma Coated Point Angel S/E



L6-084 Pt4 Titanium Plasma Coated Point STR S/E



L6-085 Pt5 Titanium Plasma Coated Point TG CVD S/E




IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان

TITANIUM PLASMA COLOUR COATED INSTRUMENTS

Hollow Handle 10mm

MIXING SPATULA WITH SHOVEL



L6-086
Bone Carrier
Titanium
185mm

L6-087
Raspator, Prichard
Length 180 mm,
10 mm hollow handle,
double-ended
10.5 mm | 4.5 mm

L6-088
Raspator, Prichard
Length 180 mm,
10 mm hollow handle,
double ended

L6-089
Raspator, Prichard
Length 180 mm,
10 mm hollow handle,
double-ended
10.5 mm | 4.5 mm

L6-090
Mixing Spatula
with shovel
170mm

L6-091
Mixing Spatula
with shovel
185mm



- Инструменты для имплантологии
- Implantologie
- Implantologie


- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



BONE GRAFTING INSTRUMENTS COMPACTOR /TAMPERS

Product offering includes instruments for harvesting, mixing and compacting bone Provides clinicians with the proper instrumentation for bone grafting procedures Weighted bone basin Prevents tipping and spilling of bone and other materials during mixing

Hollow Handle 9.5mm Length 170mm



L6-092
Graft Packer,
Concave Tips
DE 4/6mm

L6-093
Bone Spoon/
Packer DE 3mm

L6-094
Graft Packer,
Serrated Tips
DE 6/8mm

L6-095
Graft Packer,
Smooth Tips
DE 3/4mm

L6-096
Graft Packer,
Smooth Tips
DE 5/8mm

L6-097
Titanium graft packer/
placement site depth
gauge 3.3mm & 4.0mm
Diameter tips Marked
at 8, 11, 13, 15 & 18mm

L6-098
Graft Material Spoon
Concave Tip For
Larger Grafts




IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



BONE GRAFTING INSTRUMENTS

TITANIUM PLASMA COLOUR COATED Hollow Handle Ø 10 mm



L6-099

Blade and pusher
Length 180 mm,
for straight or angled,
blade: 15 x 5 mm,
pusher: 4 mm diameter
Straight

L6-100

Blade and pusher
Length 180 mm,
for straight or angled,
blade: 15 x 5 mm,
pusher: 4 mm diameter
angular

L6-101

Spoon and pusher
Length 180 mm,
for angled spoon:
10 x 5 mm, pusher:
3 mm

L6-102

Tunnelierungs
Instrument
length 175 mm,
diameter 3 mm

L6-103

Tunnelierungs
Instrument
length 175 mm,
different designs

L6-104

Tamper
Length 170 mm,
for bone material,
double-ended,
diameter 3.2 mm
(2-4-6-8-10) and diameter
2.5mm (2-4-6-8-10)



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان

BONE GRAFTING INSTRUMENTS



L6-105 Ø 4mm
L6-106 Ø 5mm
Compactor
150mm

L6-107
Sinus Window
Bone Tamper
4 mm

L6-108
Sinus Window
Bone Tamper
7 mm

L6-109
6/8mm

L6-110
3/4mm

Pluggers

Hollow Handle 9.5mm Length 160 mm, with titanium tips for mixing and pouring of autologous, heterologous bone substitute material, diameter (3/4 or 6/8 mm)






IMPLANTOLOGY

LEXER
International

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Імплантологія
- زراعة الأسنان



L6-111
10x5mm
Spoon and
Pusher
with safety
stop

L6-112
10x5mm
Spoon and
Pusher
with safety
stop

L6-113
10x5mm
Spoon and
Pusher
with safety
stop

L6-114
15x10mm


L6-115
15x10mm
Spoon and
Pusher
with safety
stop

L6-116
10x5mm
Spoon and
Pusher
with safety
stop

L6-117
10x5mm
Spoon and
Pusher
with safety
stop

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان

PERIOSTEALS ELEVATOR 10mm Hollow Handle



L6-118
Raspator,
Bone Scraper
Length 165 mm,
10 mm hollow
handle,
double-ended 6
mm | 4 mm,

L6-119
Raspator,
Raspatory
Length 180 mm,
10 mm hollow
handle,
double-ended
3 mm | 3 mm,

L6-120
Raspator,
Raspatory
Length 180 mm,
10 mm hollow
handle,
double-ended
3 mm | 3 mm,

L6-121
Raspator,
Papilla
Length 185 mm,
10 mm hollow
handle,
double-ended
3 mm | 2 mm,

L6-122
Raspator,
Bone Scraper
Length 175 mm,
10 mm hollow
handle,
double-ended,
3 mm | 2 mm,

L6-123
Raspator,
Micro-Freer
Length 180 mm,
10 mm hollow
double-ended
3 mm | 4 mm

L6-124
Raspator,
Freer
Length 180 mm,
10 mm hollow
handle,
double-ended
5 mm | 5 mm,



IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان

PERIOSTEALS ELEVATOR 10mm Hollow Handle



L6-125
Raspator,
Freer
10 mm hollow
handle,
high-quality
medical
double ended

L6-126
Raspator,
Molt
Length 180 mm,
10 mm Hollow
grip,
double-ended
6 mm

L6-127
Raspator,
Micro Freer
Length 185 mm,
10 mm hollow
handle,
double-ended

L6-128
Raspator,
Micro Freer
Length 185 mm,
10 mm hollow
handle,
double-ended

L6-129
Raspator, Micro
Raspator /
Chisel
Length 170 mm,
10 mm hollow
handle,
double-ended

L6-130
Raspator, Micro
Raspator /
Chisel
Length 175 mm,
10 mm hollow
handle,
double-ended

L6-131
Raspator,
Micro Raspator /
Chisel
Length 175 mm,
10 mm hollow handle,
double-ended
3.5mm | 3.5 mm

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان

PERIosteals ELEVATOR 10mm Hollow Handle



L6-132
Raspatory Micro Molt
Length 175 mm,
corrosion-resistant

L6-133
Raspatory Chipping
Workers:
Chisel 24G and 7.5 mm,
length about 175 mm,
corrosion-resistant

L6-134
Raspatory Molt
Work ends:
5.0 mm and 4.0 mm,
length about 175 mm,
corrosion-resistant

L6-135
Rasp, Molt M9
Length 175 mm,
corrosion-resistant

L6-136
Raspator, Wautier
Length 175 mm,
corrosion-resistant

L6-137
Raspator, Benque
Length 170 mm,
workers
5 mm | 7 mm,
corrosion-resistant

L6-138
Raspator, Nordberg
Length 180 mm,
10 mm hollow handle,
double-ended
6 mm | 4 mm



IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان

Curettes

Hollow Handle 10mm Length 170mm



L6-139 Miller



L6-140 Miller



L6-141 LUCAS #84



L6-142 LUCAS #85



L6-143 LUCAS #86



L6-144 LUCAS #87



L6-145 LUCAS #88



L6-146 LUCAS86S
SERRATED



L6-147 LUCAS86AS
SERRATED



- Инструменты для имплантологии
- Implantologie
- Implantologie

- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



If an implant cannot be set, since the diameter of the jaw is not sufficient, a piece of bone has to be implanted. The additional bone is extracted e.g. from the lower jaw. The bone piece is then locked onto the compound area with the straight forceps (front tooth sector).



TISSUE MEMBRANE FORCEPS



L6-148
Tissue Graft Forcep
Serrated Teeth
150mm

L6-149
Self Locking
Grooved
150mm

L6-150
Tissue Graft Forcep
Serrated Teeth
150mm

L6-151
Titanium Metal
Locking Plier

L6-152
TC Gripper Carbide
Tipped Tissue Forceps
Serrated Carbide Tips
145mm

L6-153
Adson Carbide
Tipped Tissue Forceps
Serrated Carbide Tips
120mm

L6-154
Chalazion Mucocele
Forceps
90mm Length
12mm x 23mm
Inside Ring

1x2



L6-155
Micor Adson
Tissue Forcep 120mm
1x2 Teeth -
0.5mm Tips

L6-156
Adson 7x7
Tissue Forceps
7x7 Teeth
120mm

L6-157
Adson 7x7 Curved
Tissue Forceps
7x7 Teeth 120mm

L6-158
Suture Pliers
Suture Forceps
with
1.6mm Ring Tip
160mm

L6-159
Membrane Forcep
With Delicate
Ring Tip
10 cm

L6-160
Debakey Atraumatic
Tissue Forceps
Straight
Serrated Teeth
160mm

L6-161
Debakey Atraumatic
Tissue Forceps Curved
Serrated Teeth
160mm



IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان

TISSUE MEMBRANE FORCEPS



L6-162
Russian Long
Tissue Forceps
Serrated Tip
200mm

L6-163
Russian
Tissue Forceps
Serrated Tip
150mm

L6-164
Tissue Punch Forceps
To Remove Tissue
After Using Punch
150mm

L6-165
Narrow
135mm

L6-166
105mm
Membrane Tweezers Kuehne

L6-167
wide
135mm

L6-168
Bishop Harman
Micro Tissue Forcep
8.6cm 1x2 Teeth
- 0.3mm Tips



L6-169
Pin-Tweezers
length 150mm
with lock

L6-170
without lock
150mm

L6-171
Haltepinzette with lock
155mm, for pins, towards
the tips, etc

L6-172
Carry Forceps
150mm


L6-173
Carry Forceps
150mm

L6-174
Carry Forceps
150mm



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان


OSTEOTOME CONVEXES / CONCAVE / CONIQUES



L6-175
SET OF 5 PIECES
Convexes

L6-176
SET OF 5 PIECES
Concaves

L6-177
SET OF 5 PIECES
Coniques



	CONVEXES STRAIGHT	L6-175	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONCAVES STRAIGHT	L6-176	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONIQUES STRAIGHT	L6-177	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm



IMPLANTOLOGY

LEXER
International

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان


OSTEOTOME CONVEXES / CONCAVE / CONIQUES



L6-178
SET OF 5 PIECES
Convexes Curved

L6-179
SET OF 5 PIECES
Concaves Curved

L6-180
SET OF 5 PIECES
Coniques Curved



	CONVEXES CURVED	L6-178	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONCAVES CURVED	L6-179	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONIQUES CURVED	L6-180	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان






IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان

BONE INJECTORS, COLLECTORS & SUCTION TUBE



L6-200
 \varnothing 2,5mm
 Straight / Bend
 180mm
 I.D 2,5mm
 O.D 3,5

L6-201
 \varnothing 3,5mm
 Straight / Bend
 180mm
 I.D 3,5mm
 O.D 4,5

L6-202
 \varnothing 4,5mm
 Straight / Bend
 180mm
 I.D 4,5mm
 O.D 5,5

L6-203
 \varnothing 6,0mm
 Straight / Bend
 180mm
 I.D 6,0mm
 O.D 7,0

L6-204
 \varnothing 2,5mm
 170mm
 I.D 2,5mm
 O.D 3,5

L6-205
 \varnothing 3,5mm
 170mm
 I.D 3,5mm
 O.D 4,5

L6-206
 \varnothing 4,5mm
 170mm
 I.D 4,5mm
 O.D 5,5

L6-207
 \varnothing 6,0mm
 170mm
 I.D 6,0mm
 O.D 7,0



L6-208
 Amalgam Gun



L6-209
 Bone Curette
 Fig 3



L6-210
 Brush Stainless Steel

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان



BONE SCRAPER STRAIGHT & CURVED

L6-211
Bone Scraper
StraightL6-212
Bone Scraper
CurvedTissue Punch
Straight

- | | |
|-----------------|-----------------|
| L6-213
2.0mm | L6-214
4.0mm |
| L6-215
5.0mm | L6-216
6.0mm |
| L6-217
7.0mm | L6-218
8.0mm |

TISSUE PUNCH STRAIGHT & CURVED




- | | |
|-----------------|-----------------|
| L6-219
2.0mm | L6-220
4.0mm |
| L6-221
5.0mm | L6-222
6.0mm |
| L6-223
7.0mm | L6-224
8.0mm |

Tissue Punch
CurvedL6-225
Bone MillL6-226
Basins For BoneL6-227
BONE CRUSHER



IMPLANTOLOGY

- Инструменты для имплантологии
- Implantologie
- Implantoologia
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان





- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان

TREPHINE

Burs with depth markings are suited for precise and simple transplantation of cylindrical bone grafts as well as for all other bone-extraction tasks.



Ø4.0 Ø5.0 Ø6.0 Ø7.0



L6-242
Tissue Punch with pin
Set of Four Pieces



L6-243
Tissue Punch
Set of Five Pieces



L6-244
Tissue Punch
Set of Five Pieces



L6-245
Trepbine Drill
Set of Eight Pieces



L6-246
Short Implant Drills

L6-247
Long Implant Drills

L6-248
Trepbine Drill
Set of Eight Pieces



IMPLANTOLOGY

LEXER
International

- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان

KEYS FOR IMPLANT SCREW



L6-249



L6-250



L6-251




L6-252

TORQUE WRENCH & RATCHET



L6-253 Star Torque Wrench 6.35mm Universal



L6-259 Star Torque Wrench 6.35mm Universal



L6-265 Star Rachet Wrench 6.35mm Universal



L6-254 Round Torque Wrench 4.0mm Universal



L6-260 Star Torque Wrench 6.35mm Universal



L6-266 Star Rachet Wrench 4.0mm Universal



L6-255 Hex Driver With 4.0mm Driver



L6-261 Star Torque Wrench 6.35mm Universal



L6-267 Round Rachet Wrench 4.0mm Universal



L6-256 Star Torque Wrench 6.35mm Universal



L6-262 Star Torque Wrench 6.35mm Universal



L6-268 Round Rachet Wrench 3.35mm Universal



L6-257 Hex Torque Wrench 4.00mm Universal



L6-263 Star Rachet Wrench 6.35mm Universal



L6-269 Star Rachet Wrench 4.0mm Universal



L6-258 Star Torque Wrench 6.35mm Universal



L6-264 Star Rachet Wrench 6.35mm Universal



L6-270 Star Rachet Wrench 6.35mm Universal



- Инструменты для имплантологии
- Implantologie
- Implantologie
- Інструменти для імплантології
- Implantoloia
- زراعة الأسنان


BONE EXPANDER SCREW & SURGICAL KITS

**L6-271**

Bone Universal Kit
 Drill: 2mm, 2.8mm, 3.2mm, 3.65mm,
 4.2mm, 4.5mm, 4.8mm, 5.2mm
 Trepheines Bit:
 3.0mm, 4.0mm, 5.0mm, 6.0mm,
 7.0mm, 8.0mm, 9.0mm, 10mm.
 Driver:
 2.5mm, 2.5mm Hand Pcs
 1.25mm, 1.25mm Hand Pcs
 Short Driver 2.5mm
 Long Driver 2.5mm
 1.25mm, Short Hex Driver Hand Pcs
 1.25mm, Long Hex Driver Hand Pcs
 Torque Wrench 4.0mm

**L6-272**

Bone Expander Kit
 Expander: 2.6mm, 3.0mm,
 3.4mm, 3.8mm, 4.3mm
 Rachet Driver
 Parallel Pin
 Pilot Drill
 Drill 1.4mm
 Drill Driver 4mm
 Hand Driver 4mm
 Rachet 4.0mm

**L6-273**

Bone Universal Kit
 Drill: 2mm, 2.8mm, 3.2mm, 3.65mm,
 4.2mm, 4.5mm, 4.8mm, 5.2mm
 Trepheines Bit:
 3.0mm, 4.0mm, 5.0mm, 6.0mm,
 7.0mm, 8.0mm, 9.0mm, 10mm.

Bone Compression:

2mm, 3.3mm, 2.8mm, 4.2mm, 3mm, 3.5mm
 Driver:
 2.5mm, 2.5mm Hand Pcs
 1.25mm, 1.25mm Hand Pcs
 Short Driver 2.5mm
 Long Driver 2.5mm
 1.25mm, Short Hex Driver Hand Pcs
 1.25mm, Long Hex Driver Hand Pcs
 Torque Wrench 6.35mm 1 Pcs

L6-274

Conical Drills: 20 Pcs
 Tissue Punches:
 3.0mm, 3.50mm, 4.0mm, 4.5mm,
 5.0mm, 6.0mm, 6.50mm
 Bone Expander:
 1.6/3.8mm, 1.9/3.0mm, 2.1/3.6mm,
 2.6/3.8mm, 2.9/4.3mm
 Drill Adapter: 1 Pcs
 1.25mm, 2.50mm Hex Driver Hand Pcs

Pilot Drill: 1.6mm

Drill Bit: 1.4mm, 1.6mm, 1.8mm, 2.0mm, 2.2mm
 Rachet Driver: 4mm/6.35, 1.25/2.50
 Rimer: 3.2/4.7mm, 3.7/4.7mm
 Parallel Pin: 4.5mm, 2.5mm
 Double step Drill: 2.4/2.8mm
 Depth Gauge: 1 Pcs
 Torque Wrench 6.35mm Universal: 1 Pcs



IMPRESSION TRAYS & HOLLOWWARE

LEXER
International

A collection of dental impression trays and hollowware, including various sizes and designs of tray sets and a tray organizer.

A smiling patient in a dental chair, ready for treatment.

A dental professional in a white coat and gloves performing a dental examination on a patient.


PREMIUM QUALITY DENTAL INSTRUMENTS

IMPRESSION TRAYS

- Оттисковые ложки
- Porte-empreinte
- Cubetas para impresiones
- صوانی الاتطباع

PERFORATED WITH RIM LOCK

Set of 06 Pcs.

**L7-001 SET OF 06 PCS**

Upper: S, M, L
Lower: S, M, L

PERFORATED WITH RIM LOCK

Set of 8 Pcs.

**L7-002 SET OF 08 PCS**

Upper: XS, S, M, L
Lower: XS, S, M, L

L7-003 SET OF 08 PCS

Upper: S, M, L, XL
Lower: S, M, L, XL

L7-004 SET OF 10 PCS

Upper: XS, S, M, L, XL
Lower: XS, S, M, L, XL



IMPRESSION TRAYS

- Оттисковые ложки
- Porte-empreinte
- Cubetas para impresiones
- صوانی الاتطبع

SOLID WITH RIM LOCK

Set of 6 Pcs



L



M



S



L



M




S

L7-005 SET OF 06 PCS

Upper: S, M, L
Lower: S, M, L

SOLID WITH RIM LOCK

Set of 8 Pcs.



L7-006 SET OF 08 PCS

Upper: XS, S, M, L
Lower: XS, S, M, L

L7-007 SET OF 08 PCS

Upper: S, M, L, XL
Lower: S, M, L, XL

L7-008 SET OF 10 PCS

Upper: XS, S, M, L, XL
Lower: XS, S, M, L, XL

IMPRESSION TRAYS

- Оттисковые ложки
- Porte-empreinte
- Cubetas para impresiones
- صوانی الاتبع

PERFORATED FULL DENTURE

Set of 8 Pcs

**L7-009 SET OF 08 PCS**


Upper: S, M, L, XL
Lower: S, M, L, XL

L7-010 SET OF 08 PCS

Upper: S, M, L
Lower: S, M, L

SOLID FULL DENTURE

Set of 8 Pcs

**L7-011 SET OF 08 PCS**

Upper: S, M, L, XL
Lower: S, M, L, XL

L7-012 SET OF 06 PCS


Upper: S, M, L
Lower: S, M, L




IMPRESSION TRAYS

- Оттисковые ложки
- Porte-empreinte
- Cubetas para impresiones
- صوانی الاتطبع

Perforated Without Rim Lock



Solid With Rim Lock



IMPRESSION TRAYS

- Оттискные ложки
- Porte-empreinte
- Cubetas para impresiones
- صوانی الاتطبع

PARTIAL SOLID



L7-049 44mm
L7-050 55mm
L7-051 60mm
L7-052 65mm



L7-053



L7-054 44mm
L7-055 55mm
L7-056 60mm
L7-057 65mm



L7-058



L7-059

Set of 11 Pcs L7-060

PARTIAL PERFORATED



L7-061 44mm
L7-062 55mm
L7-063 60mm
L7-064 65mm



L7-065



L7-066 44mm
L7-067 55mm
L7-068 60mm
L7-069 65mm



L7-070



L7-071



L7-072




L7-073

Set of 13 Pcs. L7-074



IMPRESSION TRAYS

- Оттискные ложки
- Porte-empreinte
- Abdrucklöffel
- Відбиткові ложки
- Cubetas para impresiones
- صواني الانطباع




1/2L
L7-075 1/2L/S
L7-076 1/2L/M
L7-077 1/2L/L

1/2R
L7-078 1/2R/S
L7-079 1/2R/M
L7-080 1/2R/L



1/2L
L7-081 1/2L/S
L7-082 1/2L/M
L7-083 1/2L/L

1/2R
L7-084 1/2R/S
L7-085 1/2R/M
L7-086 1/2R/L



Universal
L7-087 UN/1
Front

Universal
L7-088 UN/1
Front



L7-089 1
Front
with swivel handle

L7-090 1
Front



Stolley
L7-091
ST/1

Partial occlusal stainless steel impression trays

They allow you to take an impression
of two semi-arches antagonistic



L7-092
OCP



L7-093
OCG



COTTON HOLDER

Cotton Holder -
L7-094

L7-095

ENDO BOX



L7-096

SPIRIT LAMP



L7-097

ENDO BOX PREMIUM



L7-098

Spirit Lamp
L7-099

ENDO BOX ROUND



L7-100

L7-101
Amalgam Cup

STERILIZATION DRUM



L7-103

INSTRUMENTS TRAY WITH LID



L7-102



Made of 18/8 stainless steel


KIDNEY TRAY WITH COVER
 (English Pattern)

	Size
L7-104	150 mm
L7-105	200 mm
L7-106	240 mm
L7-107	300 mm



L7-108

**Kidney Trays**

With lid available	Without lid available
L7-109 6"	L7-113 6"
L7-110 8"	L7-114 8"
L7-111 10"	L7-115 10"
L7-112 12"	L7-116 12"



Deep
L7-117 2.25" - 2"
L7-118 4.75" - 2"
L7-119 5.5" - 2.5"
L7-120 6.25" - 2.75"
L7-121 7" - 3.5"

**Sterilizing Drum**

Art #	Ø mm	h mm
L7-122	80	80
L7-123	120	120
L7-124	160	120
L7-125	190	140
L7-126	190	200
L7-127	240	120
L7-128	240	160
L7-129	240	240
L7-130	290	120
L7-131	290	160
L7-132	290	240
L7-133	340	160
L7-134	340	240
L7-135	390	145
L7-136	390	200
L7-137	390	290



L7-138
 Endodontic Box

Forceps Jar
 L7-139 25 x 80mm
 L7-140 30 x 100mm
 L7-141 100 x 200mm
 L7-142 75 x 173mm


 L7-143 6" x 8"
 L7-144 8" x 10"
 L7-145 10" x 12"
 Without Lid

L7-146
 Instruments Box

SCALER TRAY
 L7-147
 100 x 200 mm



CASSETTE FOR
4 INSTRUMENTS
L7-148



CASSETTE FOR
5 INSTRUMENTS
L7-149



CASSETTE FOR
6 INSTRUMENTS
L7-150



CASSETTE FOR
7 INSTRUMENTS
L7-151




CASSETTE FOR
10 INSTRUMENTS
L7-152



CASSETTE FOR
20 INSTRUMENTS
L7-153



CASSETTE FOR
EXTRACTION FORCEPS
L7-154



CASSETTE FOR
ROOT ELEVATOR
L7-155



CASSETTE FOR
6 INSTRUMENTS UTILITY
L7-156



INDEX

ABOU-RASS EXPLORERS45
ACUTE IMPLANTOLOGY144
ADAMS PLIER129, 132
ADERER PLIERS132
ADHESIVE REMOVING PLIER133
ADSON 7X7 TISSUE FORCEPS160
ADSON BROWN TWEEZERS43
ADSON CARBIDE TIPPED TISSUE FORCEPS160
ADSON NEEDLE HOLDERS106
ADSON TWEEZERS43
ADSON-TC TWEEZERS43
ALBRECHT FORCEPS10
ALGINATE SPATULAS PLASTER SPATULAS & WAX KNIVES90
ALGINATE, PLASTER SPATULAS, WAX KNIVES PLASTER FOR KNIVES & PLASTER SHEAR90
ALLEN PERIOSTEALS ELEVATORS120
ALLIS 6:7 HAEMOSTATIC FORCEPS109
ALLIS-BABY 3:4 HAEMOSTATIC FORCEPS109
AMALGAM CARRIERS118
AMALGAM CUP177
AMALGAM GUN165
AMALGAM INSTRUMENTS118-119
AMALGAM WELL117, 119
AMLAGAM CONDENSER67
ANDREWS FILLING INSTRUMENTS55
ANISWORTH111
ANTERIOR59
ANTERIOR BAND REMOVER134
ANTERIORLUKS50
APEXO ERGO DYNAMIC ROOT ELEVATOR26
APEXO ROOT ELEVATOR17-19, .26, .28
APPUNTITO RETRACTORS98
ARROW CLASP BENDING PLIERS133
ARTICULATORS89
AUSTIN CHEEK RETRACTORS95
AUTOMATIC CROWN REMOVER116
BABY TOOTH EXTRACTION FORCEPS11
BACK ACTION FILLING INSTRUMENTS57
BACKHAUS TOWEL CLAMPS109
BACOLITE HANDLE90
BALDWIN FILLING INSTRUMENTS54
BALL & SOCKET PLIER133
BALL BURNISHER63
BALL FILLING INSTRUMENTS57
BAND CRIMPING PLIER131
BAND FORMING139
BAND PINCHING134
BAND PUSHER114, 136
BAND SEATER AUTO CLEAVABLE136
BAND SEATING & REMOVING INSTRUMENTS138
BAND/CROWN CRIMPING134
BASINS FOR BONE166
BAUMGARTNER NEEDLE HOLDERS107
BAUMGARTNER NEEDLE HOLDERS104
BEALE WAX AND MODELLING CARVERS88
BEAVERTAILA FILLING INSTRUMENTS57
BEEBEE CROWN SCISSORS104
BEERENDONK MINI GAUGE136
BERLIN FILLING INSTRUMENTS80
BIRD BEAK PLIERS132
BIRTLES SPATULA APPLICATOR63
BISHOP HARMAN MICRO TISSUE FORCEP161
BLADE AND PUSHER153
BLUMENTHAL BONE RONGEURS126
BLUMENTHAL MICRO BONE RONGEURS126
BOGHOSIAN COMPOSITE INSTRUMENTS77
BOGHOSIAN FILLING INSTRUMENTS54
BOLEY GAUGE136
BONE CARRIER TITANIUM151
BONE CHISEL123
BONE CRUSHER166-167
BONE CRUSHER BONE GRINDING GRAIN167
BONE CRUSHER TITANIUM METAL CUTTER167
BONE CURETTE165
BONE EXPANDER KIT170
BONE EXPANDER SCREW & SURGICAL KITS170
BONE FILE121
BONE GRAFT/CARRIER TITANIUM155
BONE GRAFTING INSTRUMENTS153-154
BONE GRAFTING INSTRUMENTS COMPACTOR /TAMPERS152
BONE INJECTORS, COLLECTORS & SUCTION TUBE165
BONE MILL166
BONE RONGEURS126
BONE SCRAPER166
BONE SPLITTING CHISELS125
BONE SPOON / PACKER152
BONE THICKNESS GAUGE141
BONE UNIVERSAL KIT170
BOON'S GAUGE137
BRACKET HOLDER TWEEZERS138
BRACKET REMOVING PLIER ANGLED TC - INSERTS129
BREWER111
BRUNINGS CHEEK RETRACTOR97
BRUSH STAINLESS STEEL165
BUCAL TUBE HOUXER138
BUCCAL DENTAL PHOTOGRAPHIC MIRRORS92
BURNISHER FILLING INSTRUMENTS54
BURNISHERS58, 78
BURNISHERS FILLING INSTRUMENTS57
BUSER PERIOSTEALS ELEVATORS120
CALIPER / CALIPERS RULE136
CANINE CONTOURING PLIER133
CAP REMOVING PLIER128
CARRY FORCEPS161
CARVERS60, 65
CARVERS INSTRUMENTS60, 64, 75
CASSETTE FOR 10 INSTRUMENTS179
CASSETTE FOR 20 INSTRUMENTS179
CASSETTE FOR 4 INSTRUMENTS179
CASSETTE FOR 5 INSTRUMENTS179
CASSETTE FOR 6 INSTRUMENTS179
CASSETTE FOR 7 INSTRUMENTS179
CASSETTE FOR EXTRACTION FORCEPS179
CASSETTE FOR ROOT ELEVATOR179
CASSETTE TRAYS179
CASTROVIEJO NEEDLE HOLDERS104, 108
CASTROVIEJO SCISSORS105
CASTROVIEJO TC NEEDLE HOLDERS108
CEMENT SPATULAS87
CERI SAW91
CHALAZION MUCOCELE FORCEPS160
CHANDLER BONE CHISEL123
CHEATEL STERILIZING FORCEPS110
CHEEK RETRACTOR97
CHEEK RETRACTORS / TISSUE RETRACTORS95
CHILDREN TOOTH FORCEPS11
CHIP BLOWER119
CHISELS73
CHOMPRET ROOT ELEVATOR30
CINCH BACK PLIER129
CLAMP, OSTEOMETER135
CLEOID-DISCOIDS FILLING INSTRUMENTS80
COGSWELL ROOT ELEVATOR22, 26
COLEMAN ROOT ELEVATOR20
COLLEGE WITH LOCK TWEEZERS42
COLLINS TONGUE FORCEPS110
COMPACTOR154
COMPASSO DE WILLIS136
COMPOSITE INSTRUMENTS64, 74, 76-79
COMPOSITE INSTRUMENTS TITANIUM BLUE74
COMPOSITE/PLASTIC FILLING INSTRUMENTS55
CONCAVES CURVED SET OF 5 PIECES163
CONCAVES SET OF 5 PIECES162
CONDENSER IMPLANTOLOGY145
CONE-SOCKET HANDLE39
CONE-SOCKET MAGNIFYING MOUTH MIRRORS40
CONE-SOCKET PLAIN MOUTH MIRRORS40
CONE-SOCKET RHODIUM39
CONICAL DRILLS: 20 PCS170
CONQUES CURVED SET OF 5 PIECES163
CONQUES SET OF 5 PIECES162
CONTOURING INSTRUMENT77
CONVEXES CURVED SET OF 5 PIECES163
CONVEXES SET OF 5 PIECES162
CORN FILLING INSTRUMENTS57
COTTON HOLDER177
COUPLAND GOUGE ROOT ELEVATOR25
COUPLAND ROOT ELEVATOR25, 29
CRANE KAPLAN TWEEZERS43
CRANE PICK ROOT ELEVATOR22
CRILE NEEDLE HOLDERS107
CRILE-WOOD NEEDLE HOLDERS106-107
CRIMPABLE HOOK PLIER TC INSERTS129
CROWN OPENER115
CROWN REMOVER114-116
CROWN REMOVER TIPS116
CROWN REMOVING PLIER115
CROWN SCISSORS & NEEDLE HOLDERS104
CROWN SPREADER114
CRUCIBLE TONGS91
CRYER ROOT ELEVATOR17, 21, 25
CURETTES84, 159
CURTIS ROOT ELEVATOR22



DARBY PERRY TRIMMER59
DE LA ROSA PLIERS132
DE LA ROSA WITH THREE GROOVES ARCH CONTOURING AND FORMING PLIER132
DEAN SCISSORS101
DEBAKEY ATRAUMATIC TISSUE FORCEPS160
DEEP BOWL178
DELICATE NEEDLE HOLDERS106-107
DELUX BOLEY GAUGE136
DENHARD CHEEK RETRACTORS94
DENTAL CAPSULE APPLIER APPLICATOR15
DENTAL PHOTOGRAPHIC MIRRORS92
DENTAL SYRINGES15
DENTAURUM RULER136
DEPTH & GINGIVAL GAUGES140
DERF NEEDLE HOLDERS104, 107
DIAGNOSTIC PROBES48
DICAL APPLICATORS FILLING INSTRUMENTS80
DIRECT BOND REMOVER (ANGLED)130
DIRECT BOND REMOVER (STRAIGHT)130
DIRECT BONDING BRACKET HOLDERS138
DISCOID CLEOID67
DISSECTING SET OF 11 PCS.127
DISSECTING SET OF 9 PCS.127
DISTAL BENDER SILICON GRIP138
DISTAL END CUTTER128
DISTAL END CUTTER FLUSH CUT128
DISTAL END CUTTER FLUSH CUT WITH SAFETY HOLD128
DISTAL IMPROVES VISIBILITY118
DOUBLE-ENDED EXPLORERS45
DRESSING FORCEPS43
DYNAMOMETER141
ELASTIC LIGATURE APPLICATOR139
ELASTIC SEPARATING PLIER111
ELLIOT MATRIX RETAINERS117
ENAMEL INSTRUMENTS72
ENDED CHEEK RETRACTORS94
ENDO BOX177
ENDO BOX PREMIUM177
ENDO BOX ROUND177
ENDO FINGER RULER SPAN MEASURE SCALE137
ENDODONTIC BOX178
ENDODONTICS INSTRUMENTS50, 62
ENDODONTICS INSTRUMENTS SOLID HANDLE51
ENDOGCLICK EXCAVATORS62
EXAMINATION KIT41
EXAMINATION KIT SET OF 3 PIECES41
EXAMINATION SET OF 4 PIECES41
EXCAVATORS62, 66, 68-69
EXCAVATORS/SURGICAL CURETTES70
EXECT BRACKET HOLDER TWEezERS138
EXPLORERS - SINGLE-ENDED44
EXPLORERS44-45
EXPLORERS INSTRUMENTS53
EXPLORERS INSTRUMENTS SOLID HANDLE45
EXPLORERS-SINGLE-END45
EXTRACTING FORCEPS (AMERICAN PATTERN)12-15
EXTRACTING FORCEPS (ENGLISH PATTERN)01-11
EXTRA-FLEX COMPOSITE INSTRUMENTS76
FARABEAU RETRACTOR95
FILLING INSTRUMENTS54-59, 63, 67, 80, 86
FILM HANGER, NAPKIN HOLDER & MOUTH PROBS127
FINGER PROTECTOR, LENGENBECK LARGE119
FINGER PROTECTOR, LENGENBECK SMALL119
FLOHR ROOT ELEVATOR17
FOERSTER DRESSING FORCEPS110
FORCEP JAR178
FORK RETRACTOR95
FREEDMAN BURNISHERS78
FREEDMAN DUCKHEAD COMPOSITE INSTRUMENTS78
FREER BONE CHISEL123
FREER PERIOSTEAL ELEVATORS122
FRET SAW91
FRIEDMAN BONE RONGEURS126
FRIEDMAN ROOT ELEVATOR25
FUNCTION OF ROOT ELEVATORS18-21
GARDNER BONE CHISEL123
GERALD TWEezERS42
GILLIES NEEDLE HOLDER WITH SCISSORS106, 108
GINGIVAL CORD PACKERS71
GOLDMAN FOX PERIOSTEAL ELEVATORS122
GOLDMAN SCALERS83
GOLDMAN-FOX PERIOSTEALS ELEVATORS120
GOLDMAN-FOX SCISSORS99, 102-103
GOLDSTEIN INTERPROXIMAL69
GRACEY CURETTES81
GRACEY SET OF 7 PCS81
GRAFT MATERIAL SPOON CONCAVE152
GRAFT PACKER152
GRAFT PACKER, CONCAVE TIPS152
GRAFT PACKER, SERRATED TIPS152
GREGG COMPOSITE INSTRUMENTS77
GREGG FILLING INSTRUMENTS67
GRITMAN WAX AND MODELLING CARVERS88
GROSS MAIER DRESSING FORCEPS110
HAEMOSTATIC FORCEPS & TOWEL CLAMPS109
HALFMoon TISSUE RETRACTOR95
HALSEY NEEDLE HOLDERS104, 107
HALSTEAD MOSQUITO 1X2 HAEMOSTATIC FORCEPS109
HALSTEAD MOSQUITO HAEMOSTATIC FORCEPS109
HALTEPINZETTE WITH LOCK161
HAMMER WITH WOODEN HANDLE91
HANDLES FOR MOUTH MIRRORS AND MIRRORS39
HARD WIRE CUTTER128, 133
HARD WIRE CUTTER 15° DEGREE ANGLE128
HARD WIRE CUTTER STANDARD STRAIGHT128
HARRIS EXTRACTING FORCEPS (AMERICAN PATTERN)13
HARTMANN (BABY MOSQUITO) HAEMOSTATIC FORCEPS109
HARTMANN RETRACTOR95
HATCHET FILLING INSTRUMENTS72-73
HATCHETS59
HEAT CARRIERS52
HEIDBRINK ROOT ELEVATOR18, 28-29
HEIDEMANN FILLING INSTRUMENTS80
HEMINGWAY CURETTES84
HEX DRIVER169
HEX TORQUE WRENCH169
HOLLENBACK CARVERS INSTRUMENTS65
HOLLENBACK COMPOSITE INSTRUMENTS76
HOLLENBACK FILLING INSTRUMENTS55
HOLLENBACK RESTORATIVE INSTRUMENTS61
HOLLOW WARE INSTRUMENTS177-178
HOOK KNIVES PLASTER SPATULAS & WAX KNIVES90
HOURIGAN PERIOSTEALS ELEVATORS120
HOURIGAN ROOT ELEVATOR25
HOW PLIER (CVD)130
HOW PLIER (STRAIGHT)130
HOWARD ROOT ELEVATOR28
HYGIENIST SCALERS83
HYGIENIST-JACQUETTE SCALERS83
HYGIENIST-TOWNER SCALERS83
HYLIN ROOT ELEVATOR20
IMPLANT DEPTH GAUGE140
IMPLANT DEPTH GAUGE & KEY140
IMPLANT GAGES140
IMPLANT SURGICAL CHISELS150
IMPLANTOLOGY140-170
IMPRESSION TRAYS171-176
IMPRESSION TRAYS PARTIAL PERFORATED SET OF 13 PCS.175
IMPRESSION TRAYS PARTIAL SOLID SET OF 11 PCS175
IMPRESSION TRAYS PERFORATED FULL DENTURE SET OF 8 PCS173
IMPRESSION TRAYS PERFORATED WITH RIM LOCK SET OF 06 PCS171
IMPRESSION TRAYS PERFORATED WITH RIM LOCK SET OF 10 PCS171
IMPRESSION TRAYS PERFORATED WITH RIM LOCK SET OF 8 PCS171
IMPRESSION TRAYS PERFORATED WITHOUT RIM LOCK174
IMPRESSION TRAYS SET OF 06 PCS173
IMPRESSION TRAYS SOLID FULL DENTURE SET OF 8 PCS173
IMPRESSION TRAYS SOLID WITH RIM LOCK174
IMPRESSION TRAYS SOLID WITH RIM LOCK SET OF 10 PCS172
IMPRESSION TRAYS SOLID WITH RIM LOCK SET OF 6 PCS172
IMPRESSION TRAYS SOLID WITH RIM LOCK SET OF 8 PCS172
INITIAL EXAMINATION SET OF 3 PIECES41
INSTRUMENTS BOX178
INSTRUMENTS TRAY WITH LID177
INTERPROXIMAL CARVER65, 77
INTERPROXIMAL CARVER, LONG77
INTERPROXIMAL OFF-ANGLE FILLING INSTRUMENTS54
INTER-PROXIMAL RESTORATIVE INSTRUMENTS61
INTRALIGAMENTAL ANESTHESIA SYRINGE15
IRIS SCISSORS99, 102-103
IVORY111
IVORY LIGHT WEIGHT111
IVORY MATRIX RETAINERS117
IWANSON SPRING CALIPER FOR METAL AND PORCELAIN135
IWANSON SPRING CALIPER FOR WAX135
JACQUETTE82
JACQUETTE SCALERS82
JARABACK PLIERS131
JOHNSON COUNTOURING PLIER129
JONES TOWEL CLAMPS109
JOSEPH SCISSORS100
KELLY HAEMOSTATIC FORCEPS109
KELLY SCISSORS100, 102-103
KELLY-RANKIN HAEMOSTATIC FORCEPS109
KEY RINGS127



INDEX

KEYS FOR IMPLANT SCREW	169	MULTIFUNCTIONAL METER CORROSION-RESISTANT	140
KIDNEY TRAY WITH COVER (ENGLISH PATTERN)	178	NAKAWALNI DENTAL LABORATORY INSTRUMENT	89
KIDNEY TRAYS	178	NANCE PLIER	131
KILNER (RAGNELL) SCISSORS	100	NAPKIN HOLDER WITH METAL CHAIN	127
KLEIN PEDODONTIC UPPER ROOT	11	NAPKIN HOLDER WITH PLASTIC CHAIN	127
KOCHER RETRACTORS	97	NEEDLE HOLDERS	104, 106, 107-108
LA GRANGE SCISSORS	100, 103	NEUMANN SCISSORS	100
LADMORE BURNISHER	63	NEVI SCALERS	82
LADMORE COMPOSITE INSTRUMENTS	78	NEVIUS EXTRACTING FORCEPS (AMERICAN PATTERN)	14
LADMORE FILLING INSTRUMENTS	57-58	NOT-PERFORATING AND NOT-DAMAGING TOWEL CLAMPS	109
LANGENBECK RETRACTORS	98	NOYES SCISSORS	105
LAXUTOR ROOT ELEVATORS SET OF 9 PCS	27	OCHSENBIN CHISELS	124
LE CRON WAX AND MODELLING CARVERS	88	OLSEN-HEGAR NEEDLE HOLDER WITH SCISSORS	106, 108
LESSMANN PLASTER SPATULAS & WAX KNIVES	90	OPTICAL	134
LIGATURE DIRECTIONAL INSTRUMENTS	139	OREGON PLUGGERS	56
LIGATURE INSTRUMENTS	139	ORINGER SELFRETAINING CHEEK-LIP RETRACTOR	93
LIGATURE TUCKER UTILITY INSTRUMENT	139	ORTHO LIGHT BAND PUSHER SCALER	138
LIGATURE TYING	134	ORTHO LIGHT DOUBLE ENDED LIGATION DIRECTOR SCALER	138
LIGATURE WIRE FORMING PLIER	129	ORTHO RULER / CONE SOCKET	141
LIGHT WIRE BIRD BEAK PLIER	132	ORTHODONTIC INSTRUMENTS	137-138
LINGUAL ARCH BENDING PLIER	131, 133	ORTHODONTICS CUTTERS TC	133
LINGUAL DENTAL PHOTOGRAPHIC MIRRORS	92	OSTEOTOME CHISELS - SET OF 4	124
LINGUAL-LARGE DENTAL PHOTOGRAPHIC MIRRORS	92	OSTEOTOME CONVEXES / CONCAVE / CONIQUES	162-163
LOCKLIN SCISSORS	101	OSTEOTOMY DEPTH GAUGE	140
LONDON COLLEGE TWEEZERS	42	P.K. THOMAS WAXING INSTRUMENTS	85-86
LORNA TOWEL CLAMPS	109	P.K.T INSTRUMENTS	85
LUCAS CURETTES	84, 159	PALMER	111
LUCAS PERIODONTAL INSTRUMENTS	124	PARMLY EXTRACTING FORCEPS (AMERICAN PATTERN)	13
LUCAS SERRATED CURETTES	159	PARTSCH BONE CHISEL	123
LUXACUT KIT	32	PARTSCH PERIODONTAL INSTRUMENTS	124
LUXATING ELEVATORS SET OF 08 PCS - BLUE TITANIUM COATED TIPS	34	PDL ELEVATORS	36, 38
LUXATING ELEVATORS SET OF 10 PCS - GOLD COATED TIPS	34	PDL ELEVATORS SET OF 10 PCS	37
LUXATING ELEVATORS SET OF 10 PCS	33	PDL ELEVATORS SET OF 7 PCS	36
LUXATING ROOT ELEVATOR	19, 26	PDL ELEVATORS SET OF 8 PCS	38
MACMILLAN ROOT ELEVATOR	25	PEDIATRIC TOOTH EXTRACTING FORCES	11
MAGNIFYING MOUTH MIRRORS	40	PERFORATEDNON-PERFORATED IMPRESSION TRAYS	176
MALLETS - PLASTER SHEARS - SAWS & BLADES	91	PERIODONTAL CHISELS	125
MALLETS	123	PERIODONTAL FILES	58
MANUAL PATTERN	115	PERIODONTAL INSTRUMENTS	120, 124-125
MARGIN INSTRUMENTS	72	PERIOSTEAL ELEVATORS - MALLETS & BONE CHISEL	123
MARKLEY PLUGGERS	56	PERIOSTEAL ELEVATORS	120, 122, 157-158, 166
MARQUETTE CONDENSER FILLING INSTRUMENTS	55	PERIOTOME GOLD TIP	35
MARQUETTE FILLING INSTRUMENTS	55	PERIOTOME LUXATING ELEVATORS SET OF 04 PCS - GOLD COATED TIPS	35
MARYLAND LASER ETCHED PROBES	46	PERIOTOMES ROOT ELEVATOR	29
MATHIEU NEEDLE HOLDERS	108	PERRY TWEEZERS	42
MATRIX RETAINERS	117	PIN LIGATION CUTTER	129
MAYO SCISSORS	101-103	PIN-TWEEZERS	161
MAYO-HEGAR NEEDLE HOLDERS	106-107	PLAIN MOUTH MIRRORS	40
MEAD BONE CHISEL	123	PLASTER KNIFE WOODEN HANDLE	90
MEAD BONE RONGEURS	126	PLASTER SHEARS	91
MEASURING INSTRUMENTS	135-136	PLASTER SPATULAS & WAX KNIVES	90
MEASURING LOOP	135	PLASTIC FILLING INSTRUMENTS	67
MEASURING LOOP, CASTROVIEJO	135	PLIERS FOR ORTHODONTICS & PROSTHETIES	128-134, 139
MEASURING PROBE CORROSION-RESISTANT	140	PLUGGERS	53, 154
MEMBRANE FORCEP	160	PLUGGERS/CONDENSERS	55-56
MEMBRANE POSITIONING IMPLANTOLOGY	145	POCKET PROBES	51
MEMBRANE TWEEZERS KUEHNE	161	POSITIVE POSITIONER	137
MERIAM TWEEZERS	42	POSTERIOR	51, 54
MERSHON SPREADER	114	POSTERIOR BAND REMOVER (LONG)	130
MESIAL/DISTAL FILLING INSTRUMENTS	58	POSTERIOR BAND REMOVER (SHORT)	130
METZENBAUM SCISSORS	101-102	POSTERIOR JACQUETTE	82
MICHIGAN PROBES	46	POTTS ROOT ELEVATOR	23
MICOR ADSON TISSUE FORCEP	160	PRICHARD PERIOSTEAL ELEVATORS	120, 122
MICRO BOLEY GAUGE STRAIGHT	141	PROBES	46-47
MICRO BOLEY GAUGE TWISTED	141	PROBES & PERIODONTAL INSTRUMENTS	49
MICRO BONE RONGEURS	126	PROXIMATOR ROOD ELEVATOR	37
MICRO CASTROVIEJO SCISSORS	105	RAGNELL - DAVIS RETRACTORS	98
MICRO FRIEDMAN MICRO BONE RONGEURS	126	RAMUS RETRACTORS	98
MICRO IRIS SCISSORS	99	RASP, MOLT	158
MICRO LIGATION CUTTER	128	RASPATOR, BENQUE	158
MIDDLEDORPF RETRACTORS	98	RASPATOR, BONE SCRAPER	166
MILLER - COLBURN BONE FILE	121	RASPATOR, FREER	157, 166
MILLER APEXO ROOT ELEVATOR	23	RASPATOR, MICRO FREER	157
MILLER BONE FILE	121	RASPATOR, MICRO RASPATOR / CHISEL	157
MILLER TWEEZERS	43	RASPATOR, MICRO-FREER	166
MILLER-COLBURN BONE FILE	121	RASPATOR, MOLT	157
MINI FRIEDMAN BONE RONGEURS	126	RASPATOR, NORDBERG	158
MINI LIGATION CUTTER	129	RASPATOR, PAPILLA	166
MINI RETRACTOR	97	RASPATOR, PRICHARD	151
MIXING SPATULA WITH SHOVEL	151	RASPATOR, RASPATORY	166
MOLT PERIOSTEAL ELEVATORS	120, 122	RASPATOR, WAUTIER	158
MORTONSON FILLING INSTRUMENTS	55	RASPATORY CHIPPING WORKERS CHISEL	158
MOUTH GAG MOLT	94	RASPATORY MICRO MOLT	158
MOUTH MIRRORS	40, 49	RASPATORY MOLT	158
MOUTH PROPS	127	READ ROOT ELEVATOR	20
MOVEABLE HEIGHT GAUGE	137	RESTORATIVE INSTRUMENTS	57, 59, 61
MULTIFUNCTIONAL METER	140	RESTORATIVE INSTRUMENTS SOLID HANDLE	56



RESTORATIVE TITANIUM INSTRUMENTS	76-79
RETRACTORS - MOUTH GAGS	94
RETRACTORS - MOUTH GAGS COLUMBIASTERNBURG	93
RETRACTORS	95, 97-98
REYNOLDS CONTOURING PLIERS	132
RHODIUM MIRROR	39
ROACH WAX AND MODELLING CARVERS	88
ROCHESTER-KOCHER 1X2 HAEMOSTATIC FORCEPS	109
ROCHESTER-PEN HAEMOSTATIC FORCEPS	109
ROGGE STERILIZING FORCEPS	110
ROMEROWSKI	58
ROOT CANAL EXPLORER	62
ROOT CANAL PLUGGER / SPREADER	50-53
ROOT ELEVATOR	17-38
ROSER-KONIG CHEEK RETRACTORS	94
ROUND RACHET WRENCH	169
ROUND TORQUE WRENCH	169
ROUNDED-CONE FILLING INSTRUMENTS	54
RUBBER DAM CLAMPS	113
RUBBER DAM CLAMPS KITS	112
RUBBER DAM FRAME	111
RUBBER DAM INSTRUMENTS	111
RUBBER DAM KITS	112
RULER	141
RULER WEDGE	136
RUSSIAN LONG TISSUE FORCEPS	161
RUSSIAN TISSUE FORCEPS	161
SCALER TRAY	178
SCALERS	82-83
SCALPEL HANDLES	93
SCHLUGER FILLING INSTRUMENTS	58
SCHMECKEBIER APEXO ROOT ELEVATOR	23
SCHWARZ PLIERS	133
SCISSORS	99, 100-103, 105
SCREENING PROBE PROBES	46
SELDIN PERIosteal ELEVATORS	122
SELDIN ROOT ELEVATOR	17-18, 22
SELF LOCKING GROOVED	160
SEMKEN TAYLOR TWEEZERS	42
SENN-MILLER RETRACTORS	98
SEPARATING MODULE	134
SHEPHERD'S HOOK EXPLORER	45
SICKLE SCALERS	83
SINGLE ENDED PROBE SOLID HANDLE	44
SINUS IMPLANTOLOGY	144
SINUS LIFT CURETTES AND ELEVATORS	148
SINUS LIFT ELEVATOR	146, 149
SINUS LIFT INSTRUMENT	145-147, 149
SINUS OBTUSE IMPLANTOLOGY	144
SINUS WINDOW BONE TAMERS	154
SINUSLIFT IMPLANTOLOGY	145
SIQVELAND RETAINERS	117
SLIM BIRD BEAK PLIER	132
SPECIAL INOX PLASTER SPATULAS & WAX KNIVES	90
SPENCER SCISSORS	101
SPRIT LAMP	177
SPoon AND PusHER	153
SPoon AND PusHER (TRINOVO)	155
SPoon AND PusHER WITH SAFETY STOP	155
SPREADERS	52
STANDARD LIGATURE CUTTER	129
STAR RACHET WRENCH	169
STAR RACHET WRENCH 6.35MM UNIVERSAL	169
STAR TORQUE WRENCH	169
STAR TORQUE WRENCH 6.35MM UNIVERSAL	169
STEP PLIER	133
STERILIZATION DRUM	177
STERILIZING FORCEPS	110
STIEGLITZ ROOT SPLINTER FORCEPS	10
STOCKES	111
STOLLEY IMPRESSION TRAYS	176
STRIP SAW	136
SUCTION CANNULA KIT	164
SUCTION NEEDLE	164
SUCTION PIPE	164
SUCTION TUBES	164
SUGARMAN FILE FILLING INSTRUMENTS	58
SUTURE FORCEPS	160
SUTURE PLIERS	160
SYNDESMOTOME ROOT ELEVATOR	37
T. BALL BURNISHER	63
TAMPER FOR BONE MATERIAL	153
TANNER COMPOSITE INSTRUMENTS	76
TANNER FILLING INSTRUMENTS	57, 67
TANNER RESTORATIVE INSTRUMENTS	61
TC GRIPPER CARBIDE TIPPED TISSUE FORCEPS	160
TEMPORARY	115
THARP CARVERS INSTRUMENTS	65
THE 2 GLICK ENDO EXACVATOR	70
THE 31L ENDO EXACVATOR	70
THE 31W ENDO EXACVATOR	70
THE 32L ENDO EXACVATOR	70
THE 33L ENDO EXACVATOR	70
THESE MULTIPURPOSE SPECIALTY PLIERS	131
THIN SPREADERS	52
THINN WILLIAMS PROBES	46
THINNER SPREADERS	52
THOMPSON CARVERS INSTRUMENTS	65
TISSUE FORCEPS	43
TISSUE GRAFT FORCEP	160
TISSUE MEMBRANE FORCEPS	160-161
TISSUE PUNCH	166
TISSUE PUNCH FORCEPS	161
TISSUE PUNCH SET OF FIVE PIECES	168
TISSUE PUNCH STRAIGHT & CURVED	166
TISSUE PUNCH WITH PIN SET OF FOUR PIECES	168
TITANIUM GRAFT PACKER	152
TITANIUM METAL LOCKING PLIER	160
TITANIUM PLASMA COLOUR COATED INSTRUMENTS	151
TOFFLEMIRE-UNIVERSAL MATRIX RETAINERS	117
TORQUE BENDING PLIER	130
TORQUE WRENCH & RATCHET	169
TOWNER SCALERS	83
TOWNER-JACQUETTE SCALERS	83
TREPILINE	168
TREPILINE DRILL SET OF EIGHT PIECES	168
TRIMMERS INSTRUMENTS	72
TRIPLEX SYRINGES SET	164
TRIPLEX SYRINGES TIPS	164
TUNNELIERUNGS INSTRUMENT	153
TURKEY PATTERN	115
TWEED ARCH FORMING PLIER	131
TWEED LOOP FORMING PLIERS	131
TWEED PLIER SLIM	131
TWEEZERS	42-43
U OF MINNESOTS CARVERS INSTRUMENTS	65
UNIVERSAL CROWN SCISSORS	104
UNIVERSAL FILLING INSTRUMENTS	54
UNIVERSITY OF NORTH CAROLINA EXPLORERS	45
UNIVERSITY OF NORTH CAROLINA PROBES	46
UNIVERSITY OF MINNESOTA RETRACTOR	95
V BEND PLIER	132
VEHE WAX AND MODELLING CARVERS	88
WAKAI SPREADERS	52
WAKEFIELD PERIODONTAL INSTRUMENTS	124
WALL RESTORATIVE INSTRUMENTS	61
WANGEN HOLDER	93
WARD PLUGGERS	56
WARD RESTORATIVE INSTRUMENTS	61
WARWICK JAMES ROOT ELEVATOR	29
WAX AND MODELLING CARVERS	88
WEIDER CHEEK & TONGUE DEPRESSOR	95
WEINGARD PLIER SLIM	130
WEINGARD PLIER UTILITY	130
WEINGART TC BAND FORMING	139
WEITLANER STAINLESS STEEL DOUBLE	94
WEST PERIosteal ELEVATORS	122
WESTCOTT SCISSORS	105
WESTON CEMENT SPATULAS	87
WILAND MINN. RESTORATIVE INSTRUMENTS	61
WILLIAM PROBE SET DOUBLE ENDED	48
WILLIAMS PROBES	46
WILLIGER PERIosteal ELEVATORS	120, 122
WINTER EXTRACTING FORCEPS (AMERICAN PATTERN)	12
WINTER ROOT ELEVATOR	31
WOODEN HANDLE PLASTER SPATULAS & WAX KNIVES	90
WOODWARD HAND PRESSURE CHISEL ROOT ELEVATOR	25
WOODWARD ROOT ELEVATOR	24
X RAY FILM HANGERS	127
YOUNG PLIER	131
YOUNG TONGUE FORCEPS	110
ZAHLE WAX AND MODELLING CARVERS	88
ZARRAL NEEDLE HOLDERS	108
ZURICHER GAUGE	136



Manufacturers, & Exporters of Dental, Surgical, Orthodontic, Hollow Wears & Manicure Instruments.

Wazirabad Road, Harrar,
Near Total Petrol Pump
Sialkot 51310 - Pakistan.

Ph : +92 52 3242644
Ph : +92 52 3242645
Fax : +92 52 3574225

E-mail: info@lexerintl.com
lexer_intl@yahoo.com
website: www.lexerintl.com



ISO 9001:2015 EN ISO 13485:2012
AC:2012